## Instructions For Use

### 1. Intended Use

The Orthodontic Bracket is a device intended for use in orthodontic treatment. The device is affixed to a tooth so that pressure can be exerted on the teeth.

### 2. Indication for use

Adhesive flash should be removed at bonding in order to minimize staining of the adhesive and to reduce the risk of bracket breakage during the debonding procedure. No special primers or pretreatments are necessary due to the bracket's special bonding base. Orthodontic brackets can be used with traditional direct or indirect bonding methods. For bonding procedures please follow the adhesive manufacturer's recommendations. When placing orthodontic brackets, it is suggested that the bracket be placed in a sliding motion, occlusal to gingival, forcing excess adhesive to the incisal edge.

### 3. Warnings

- To ensure proper function of the bracket tie-wings and avoid damage to bracket, do not strike the bracket with an instrument while attaching and removing elastomeric or metal ligature ties.
- Must dispose of and not reuse detached product. "SINGLE USE ONLY".
- Must follow orthodontist's direction if using for a different application.
- Instruct patients not to chew or bite on hard substances such as hard candy, ice, carrots, etc.
- If an allergic reaction occurs, direct patient to consult their dentist.
- It is recommended to bond these orthodontic brackets using a light curing adhesive.

### 4. Contraindications

- Heavy smoking, tobacco and alcohol abuse.
- Systemic blood disorders.
- Uncontrolled diabetes.

### 5. How to use

- Apply orthodontic adhesive on the bracket's base and attach it to the teeth.
- Insert a wire into bracket's slot.
- During treatment, patients and staff must follow the instructions of the orthodontist.

### 6. Storage condition

- Keep at room temperature.
- Keep in package or airtight box.
- Do not impact.

### 7. Debonding Instructions

- Remove flash from around the brackets with a high speed debonding bur.
- The bracket is designed to release at the bracket-adhesive interface.
- It will leave a layer of adhesive on the enamel for safety. This means that most of the adhesive will remain on the tooth after debonding.
- Remove remaining adhesive with a high speed debonding bur or adhesive removing plier.

Manufactured for: (U.S. Only)  
 Ortho Organizers, Inc.  
 1822 Aston Avenue,  
 Carlsbad, CA 92008-7306  
 USA  
 Tel: 800.547.2000, +(1) 760 448 8600  
 Fax: 800.888.7244, +(1) 760 448 8607  
 USASales@HenryScheinOrtho.com  
 IntlOrders@HenryScheinOrtho.com

mdi Europa GmbH  
 Langenhagener Str. 71  
 30855 Langenhagen, Germany  
 49 511 39 08 9530

Symbols glossary provided electronically at [HenryScheinOrtho.com/symbols](http://HenryScheinOrtho.com/symbols)

AW424-671 Rev B

# NeoCrystal Plus™

Sapphire Ceramic Bracket System

ORTHO  
 Organizers  
 A HENRY SCHEIN® COMPANY

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