EVEREST™

Orthodontic Adhesive

Light-Cure Orthodontic Adhesive

For professional use only

Instructions for Use

INDICATIONS:

Everest Light-Cure Orthodontic Adhesive is designed for use as an adhesive for retaining orthodontic brackets onto tooth enamel.

DESCRIPTION:

Everest Light-Cure Orthodontic Adhesive is an orthodontic bracket adhesive system. It includes a thick viscosity composite adhesive, a low viscosity resin sealant and an etchant gel. The light-curing mechanism assures a quick setting to a strong adhesive system.

FORMULATION:

Everest Light-Cure Orthodontic Adhesive is a composite based on a BisGMA resin system with a small particle glass filler and a photosensitive catalyst system.

The sealant is an unfilled resin system based on BisGMA with a photosensitive catalyst system.

The etch gel is a phosphoric acid gel.

DIRECTIONS FOR USE:

PREPARATION

- Prophy the teeth with a non-fluoride, oil-free prophy paste. Wash with an oil-free air/ water spray. Do not allow the patient to rinse.
- Isolate teeth with rubber dam or cotton rolls.
- Blow dry completely with an oil/moisture-free air source.

ETCHING

- Dispense 1 or 2 drops of the etch gel into a mixing well. One may also apply the etch gel directly to the tooth with a syringe and a syringe tip. Recap the etch gel.
- Apply the etch gel to the enamel with a brush or sponge pellet. Dab the etch gel onto the tooth, but do not rub against the enamel. Allow etch gel to remain for 30 seconds.
- With a few drops of water from the air/water syringe and with high volume suction near the tooth wash away the etch gel. This avoids splashing the etch. After the etch gel is gone continue to wash the area of enamel completely for 20 seconds with an air/ water spray and adequate suction. Do not allow saliva to contaminate etched enamel surfaces. Contaminated surfaces need to be re-etched for 10 seconds.
- Air dry thoroughly using oil/moisture-free air. After drying the etched surfaces of teeth should have a dull, white, frosted appearance. If the enamel is not chalky appearing, repeat the process allowing etch gel to remain on the teeth for additional 30 second increments of time and wash and dry as described in this etching procedure. This repetition may be needed especially on teeth high in fluoride content.

SEALING

- Dispense 2 drops of sealant liquid into a clean mixing well. Recap the sealant bottle.
- Apply the sealant directly onto each etched, dry surface. Remove excess resin sealant with a brush or a gentle stream of oil-free, moisture-free air. Avoid contact with gingiva or other soft tissue.

- Curing the sealant is not necessary at this time but can be done to limit the mobility of the bracket. Hold the end of the light guide as close to the tooth as possible without touching the resin. Expose the Sealant for 10 seconds on each tooth. **WARNING:** Use a blue curing light (bandwidth of 470-480 nanometers with a minimum of 500 mW/cm² intensity.)

ADHERING

- Apply a small increment of adhesive paste directly onto the bracket back, or use a small, flat ended hand instrument to place the composite on the bracket. Make sure the paste is worked into the bracket backing. Recap the syringe or syringe tip after use.

EVEREST™

Orthodontic Adhesive

Manufactured for: (US Only) Ortho Organizers, Inc. 1822 Aston Ave., Carlsbad, CA 92008 Ph: 760.448.8600/800.547.2000 HenryScheinOrtho.com

Hecho en los EEUU/ Fabriqué aux Etats-Unis
Hergestellt in den USA/ Prodotto negli USA

- Seat the bracket onto the tooth with light pressure and position the bracket. Press firmly to place but do not express all of the paste from under the bracket.
- Remove any gross excess of paste that expresses from around the bracket.
- Light-cure the adhesive from two directions (incisal and gingival or mesial and distal) for 20 seconds each direction for each bracket. The light guide tip should be within 5 mm of the tooth surface. **WARNING:** Use a blue curing light (bandwidth of 470-480 nanometers with a minimum of 500 mW/cm² intensity.)
- Arch wires can be placed immediately after bonding the last bracket.

BRACKET REMOVAL

At completion of the orthodontic treatment, remove the bracket by placing the cutter edges at the adhesive/tooth interface. Remaining adhesive can be removed with a rubber wheel and/or scaler.

WARNINGS:

- Dentists and assistants should wear gloves, masks and protective eye wear. Patients also should wear eye protection.
 - If contact with eyes occurs, flush immediately and for a prolonged period with copious amounts of water and consult a physician.
 - If contact with skin occurs, wash with soap and copious amounts of water.
 - In case of ingestion, seek medical attention immediately.
- Etch gels are caustic and can cause tissue damage in contact with skin or eyes.
- Uncured monomers in this resin material may cause an allergic reaction or skin irritation in some individuals. Do not use in individuals with allergies to these materials.
- Any non-sterilizable item that is handled in the delivery of the dental service, should be disinfected by standard dental office hygiene procedures.
- Do not allow saliva or water contamination of etched tooth surface or unset material.
- If a fluoride treatment is planned at the same appointment, apply the fluoride after bonding the bracket, and the adhesive is completely set.
- Test the dental curing light periodically to assure sufficient intensity of the light output.
- For additional information consult the Safety Data Sheet.

STORAGE:

- The system is designed to be stored between 10 °C - 24 °C (50 °F - 75 °F).
- Replace caps immediately after dispensing.
- Bring to room temperature prior to use.
- Do not store materials in the proximity of eugenol-containing products.
- Do not store the materials in intense light, or under wet conditions.

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0470 mdi Europa GmbH