### No-Mix Adhesive

**For Professional Use Only**

#### Instructions for Use

INDICATIONS:

No-Mix Adhesive is designed for use as an adhesive for retaining orthodontic brackets onto tooth enamel.

DESCRIPTION:

No-Mix Adhesive is a chemically cured orthodontic bracket adhesive system which does not require any premixing procedure. It includes a thick viscosity composite adhesive paste, a low viscosity liquid resin primer and an etch gel. The application procedure produces a quick setting, strong bracket adhesive system.

FORMULATION:

No-Mix Adhesive Paste is a composite based on a BisGMA resin system with a quartz filler and benzoyl peroxide initiator.

No-Mix Adhesive Primer is an unfilled resin system based on BisGMA with an amine catalyst.

The etch gel is a phosphoric acid gel.

DIRECTIONS FOR USE:

PREPARATION

- Prophy the teeth with a non-fluoride, oil-free prophy paste. Wash with an oil-free air/water spray. Do not allow the patient to rinse.
- Isolate teeth with rubber dam or cotton rolls.
- Blow dry completely with an oil/moisture-free air source.

ETCHING

- Dispense 1 or 2 drops of etch gel onto the mixing pad. One may also apply the etch gel directly onto the enamel with a syringe and a syringe tip. Recap the etch gel.
- Apply the etch gel to the enamel with a brush or sponge pellet. Dab the etch gel onto the tooth, but do not rub against the enamel. Allow etch gel to remain for 30 seconds.
- With a few drops of water from the air/water syringe and with high volume suction near the tooth wash away the gel. This avoids splashing the etch gel. After the etch gel is gone continue to wash the area of enamel completely for 20 seconds with an air/water spray and adequate suction. Do not allow saliva to contaminate etched enamel surfaces. Contaminated surfaces need to be re-etched for 10 seconds.
- Air dry thoroughly using oil/moisture-free air. After drying the etched surfaces of teeth should have a dull, white, frosted appearance. If the enamel is not chalky appearing, repeat the process allowing etch gel to remain on the teeth for additional 30 second increments of time and wash and dry as described in this etching procedure. This repetition may be needed especially on teeth high in fluoride content.

#### NO MIX PASTE/PRIMER APPLICATION

- Dispense a small amount of Primer into a clean mixing well. Recap the Primer bottle.
- Apply the Primer directly onto each etched, dry surface. Remove excess Primer with a brush or a gentle stream of oil-free, moisture-free air. Avoid contact with gingiva or other soft tissue.
- Apply a thin coat of Primer onto the bracket backing and thin as previously described.
- Using the preposed tip in a dispenser, apply the paste directly onto the bracket backing. If using a syringe of paste, dispense on pad and apply to the bracket backing. Recap after use.

**EVEREST™**  
Orthodontic Adhesive

Manufactured for: (US Only)  
Ortho Organizers, Inc.  
1822 Aston Ave., Carlsbad, CA 92008  
Ph: 760.448.8600/800.547.2000  
HenryScheinOrtho.com

**4070**

**EC REP**  
**mdi Europa GmbH**  
Langenhagener Str. 71  
30855 Langenhagen, Germany

- Seat and position the bracket onto the tooth within 30 seconds, pressing firmly to place but do not express all of the paste from under the bracket. Do not attempt to reposition brackets after this time.
- Remove any gross excess of paste that expresses from around the bracket.
- Arch wires can be placed 5 minutes after bonding the last bracket.

BRACKET REMOVAL

- At completion of the orthodontic treatment remove the bracket by placing the cutter edges at the adhesive/tooth interface. Remaining adhesive can be removed with a rubber wheel and/or scaler.

WARNINGS:

- Dentists and assistants should wear gloves, masks and protective eye wear. Patients also should wear eye protection.
  - If contact with eyes occurs, flush immediately and for a prolonged period with copious amounts of water and consult a physician.
  - If contact with skin occurs, wash with soap and copious amounts of water.
  - In case of ingestion, seek medical attention immediately.
- Etch gels are caustic and can cause tissue damage in contact with skin or eyes.
- Uncured monomers in this resin material may cause an allergic reaction or skin irritation in some individuals. Do not use in individuals with allergies to these materials.
- Any non-sterilizable item that is handled in the delivery of the dental service, should be disinfected by standard dental office hygiene procedures.
- Do not allow saliva or water contamination of etched tooth surface or unset material.
- If a fluoride treatment is planned at the same appointment, apply the fluoride after bonding the bracket, and the adhesive is completely set.
- For additional information consult the Safety Data Sheet.

STORAGE:

- The system is designed to be stored between 35 °F - 50 °F (2 °C - 10° C).
- Bring to room temperature prior to use.
- Replace caps immediately after dispensing.
- Do not store materials in the proximity of eugenol-containing products.
- Do not store the materials in intense light, or under wet conditions.

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