## Glass Ionomer Luting Cement

For professional use only

Instructions for Use

### INDICATIONS:

Glass Ionomer Luting Cement is used for the cementation of orthodontic bands. It is also an effective luting agent for conventional crowns and bridges. Other uses may include base or lining under all restorative materials or following root canal treatment. It may also be used as a temporary filling material.

### DESCRIPTION:

Glass Ionomer Luting Cement is a naturally fluoride releasing cement that may reduce the incidence of decalcification and caries. Glass Ionomer Luting Cement chemically bonds to enamel and dentin without the use of any etchants.

### FORMULATION:

Powder: fluoro-alumina-silica glass

Liquid: water solution of polyacrylic acids

### DIRECTIONS FOR USE:

1. Prophyl the teeth with flour of pumice. Do not use a prepared prophyl paste.
2. Wash and isolate the teeth to be banded. Blow the excess water off of the teeth but do not dehydrate the tooth surface.
3. Remove the top sheet of the mixing pad to provide a fresh clean sheet for mixing on the remaining pad. One may also use a glass slab either cooled or frozen to extend working time. A frozen slab must be wiped clean of any water condensation.
4. Shake the powder bottle and tap it on a firm surface to properly disperse and condense the powder.
5. Measure out scoops of powder with the scoop provided. Wipe off excess powder with a spatula to make the powder flush with the top of the scoop. Eject powder onto mixing pad or glass slab by turning upside down and tapping the side of the measuring scoop with a spatula.
6. Dispense liquid drops. Always hold the bottle completely vertical when dispensing drops. Gently squeeze bottle to produce drop and release pressure after each drop is dispensed. The ratio of powder/liquid depends on the bottle in the kit as follows:
  - a. For 50 g powder / 50 ml liquid kit: 1 scoop of powder to 1 drop of liquid;
  - b. For 15 g powder / 15 ml liquid kit: 2 scoops of powder to 3 drops of liquid;
  - c. For testing purposes: 1.6 g of powder to 1 g liquid.
7. Divide the powder into three equal amounts and mix the first increment for 10 seconds, then incorporate the second for 10 seconds, followed by the last increment for 10 seconds. Mix quickly until a creamy consistency is reached.
  - a. Mixing time -- 30 seconds
  - b. Working time -- 1 minute
8. Apply the mixed cement to the inside of the band immediately after completing the mix and seat the band while the cement is still glossy. Keep the teeth isolated until the cement gels.
9. Excess material should be removed before the cement is completely hardened.
  - a. Final setting time -- 6 to 7 minutes.

### BAND AND CEMENT REMOVAL:

Remove the bands with any pliers designed for that purpose. After the band has been removed, the excess cement left on the tooth can be removed by drying with a direct stream of air. This desiccation of the cement makes the cement easier to remove.

### WARNINGS:

- Dentists and assistants should wear gloves, masks and protective eye wear. Patients also should wear eye protection.
  - If contact with eyes occurs, flush immediately and for a prolonged period with copious amounts of water and consult a physician.
  - If contact with skin occurs, wash with soap and copious amounts of water.
  - In case of ingestion, seek medical attention immediately.
- Do not use in individuals with allergies to these materials.
- Any non-sterilizable item that is handled in the delivery of the dental service, should be disinfected by standard dental office hygiene procedures.
- Do not allow saliva or water contamination of material during the placement procedure.
- Due to the chemically adhesive nature of Glass Ionomer Luting Cement, all reusable instruments that come into contact with the cement, i.e., spatulas, glass slabs, and placement instruments, should be wiped clean under running water immediately after use.
- For additional information consult the Safety Data Sheet.

### STORAGE:

- The system is designed to be stored between 10 °C - 24 °C (50 °F - 75 °F). Do not refrigerate.
- Replace caps immediately after dispensing.

Manufactured for: (US Only)  
Ortho Organizers, Inc.  
1822 Aston Ave., Carlsbad, CA 92008  
Ph: 760.448.8600/800.547.2000  
HenryScheinOrtho.com

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**mdi Europa GmbH**