## Light-Cure Orthodontic Band Cement

[English]

For professional use only

Instructions for Use

### INDICATIONS FOR USE:

Band Cement is intended for the cementation of orthodontic bands onto posterior teeth.

### DESCRIPTION:

Band Cement is a one paste, light-curing, fluoride releasing, glass ionomer/resin cement.

### FORMULATION:

Band Cement contains a Bis-GMA based resin system with a chemical catalyst system.

### DIRECTIONS FOR USE:

Bring to room temperature prior to use.

1. Roughen the inside of the band with a fine diamond bur or air abradar to improve the adhesion of the cement to the band, if the band is not micro-etched by the manufacturer.
2. Prepare teeth with a non-fluoride, oil-free prophyl paste. Wash with an oil-free air/water spray.
3. Isolate teeth with cotton rolls or rubber dam. **Nota:** Etching enamel is not necessary but may increase retention with high stress appliances.
4. Blow dry completely with an oil/moisture free air source.
5. Dispense Band Cement onto a mixing pad or directly on the band. It is not necessary to fill the band with cement. **Warning:** Protect the dispensed paste from ambient light. Do not place the material under the operator light. Exposure to these lights may cause a premature setting of the cement.
6. Replace cap on the syringe immediately after dispensing.
7. Seat the band(s).
8. Clean off excess cement. A one-second exposure with the curing light can cause a rubbery set that can facilitate the removal of the excess cement.
9. Light cure for 30 seconds from occlusal direction using a light with a minimum output of 500 mW/cm<sup>2</sup>. With a high intensity light (approximately 2 W/cm<sup>2</sup>) exposures of 6 seconds will cure the material. Curing for this time from both occlusal and gingival directions will strengthen the set of the material. Cure the cement thoroughly with the light tip close to the material. The band can be loaded immediately after light curing.

### Other Applications:

- a) Band Cement can be used to bond brackets, lingual retainers, etc. When bonding these appliances, you first must etch the enamel, and a light curing sealant should be applied to the enamel.
- b) Band Cement can retain to porcelain surfaces with the use of a porcelain etch and porcelain primer. The porcelain may require a slight roughening with a diamond bur to enhance the etching process.
- c) Band Cement can also bond to plastic with the aid of a plastic primer.

### WARNINGS:

- Dentists and assistants should wear gloves, masks and protective eyewear. Patients also should wear eye protection.
  - If contact with eyes occurs, flush immediately and for a prolonged period with copious amounts of water and consult a physician.
  - If contact with skin occurs, wash with soap and copious amounts of water.
  - In case of ingestion, seek medical attention immediately.
- Uncured monomers in this resin material may cause an allergic reaction or skin irritation to some individuals. Do not use in individuals with allergies to these materials.
- Any non-sterilizable item that is handled in the delivery of the dental service should be disinfected by standard dental office hygiene procedures.
- Do not allow saliva or water contamination of tooth surface or unset material.
- Test the dental curing light periodically to assure sufficient intensity of the light output.
- For additional information consult the Safety Data Sheet.

### STORAGE:

- The system is designed to be stored between 59 °F - 75 °F (15 °C - 24 °C).
- Do not store materials in the proximity of eugenol-containing products.
- Do not expose materials to intense light or elevated temperatures.

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Ortho Organizers, Inc.
**mdi Europa GmbH**