[English]

## Dispenser for Tips

*For Professional Use Only*

*Instructions for Use*

### INDICATION:

The Dispenser for Tips is designed to be used to deliver dental composite resin restoratives and cements from unit dose tips to tooth preparations.

### DESCRIPTION:

- The dispenser is designed to express composite material easily from a unit dose tip.
- The dispenser provides easy control of the amount of material expressed from the tip.
- The dispenser has a loading end that can accept most of the unit dose tips in the dental market.
- The dispenser allows a more ergonomic delivery of the material into the desired location.
- The dispenser allows a more open visual field into the preparation of the tooth.
- The dispenser is autoclavable to allow proper asepsis.

### DIRECTIONS FOR USE:

1. Sterilize before use. Follow normal sterilization procedures for steam heat autoclaving, chemclave, or dry heat oven.
2. Insert a loaded tip into the opening of the dispenser.
3. Snap the tip to place into the dispenser.
4. Rotate the tip to the desired position.
5. Squeeze the handle slowly to exude the proper amount of material. The force needed will vary with the viscosity of the material.

### WARNINGS:

- Do not exceed 165 °C when sterilizing.
- Do not use when hot from the sterilization procedure.
- Examine the dispenser for suitability before use.

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### Dispenser per Puntali

 **Manufacturer**  
**Made in the USA**

Manufactured for: (US Only)  
Ortho Organizers, Inc.  
1822 Aston Ave., Carlsbad, CA 92008  
760.448.8600/800.547.2000  
HenryScheinOrtho.com

 **EC REP**  
**mdi Europa GmbH**  
Langenhagener Str. 71  
30855 Langenhagen, Germany

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