INSTRUCTIONS FOR USE:

1. Position the CFI Opening Tool with the Double-V-Pronged end toward the gingival and the Single Point toward the occlusal (Fig. 1).
2. Center the CFI Opening Tool perpendicularly to the face of the bracket. Place the Single Point into the funnel located on the face of the SLX Bracket. Then seat the Double-V-Pronged end onto the archwire. Keep a loose grip on the tool, as activation of the hinge will alter the fit (Fig. 2).
3. While supporting the patient's tooth (Fig. 1), squeeze the tool fully. This will open the bracket door (Fig. 3). Immediately remove the tool from the bracket.


Note: To close the Carriere SLX Bracket, use light finger pressure to slide the door gingivally. Listen, or feel, for one click as the door starts to close. A second click will be heard or felt as the external clasp and spring fully close over the Sure-Lok™ Clip (Fig. 4).


(Fig. 1)

(Fig. 2)

(Fig. 3)

(Fig. 4)

REPROCESSING ORTHODONTIC INSTRUMENTS

- Dental instruments—for intended use only by dental professionals.
- Warning: This product contains nickel and chromium and should not be used for individuals with known allergic sensitivity to these metals.
- All instruments must be inspected, cleaned, and sterilized before each use.
- The life time of instruments depends on the frequency of use, the care of the user, and proper reprocessing methods.
- The user is responsible for inspecting instruments prior to each use, and for the use of damaged and dirty instruments (Ortho Organizers, Inc. has no liability in case of disregard).

Cleaning

- Water quality may influence the result of the cleaning and disinfection of the instruments. Corrosion could be caused by high contents of chloride or other minerals in tap water. If problems with stains and corrosion occur and all other reasons can be excluded, it might be necessary to test the tap water quality in your area. With the use of completely deionized or distilled water most water quality problems can be avoided beforehand.
- Clean and remove any bulk debris immediately after use. Automated washing is not recommended. Do not clean any instruments using metal brushes or steel wool. It is recommended to use an enzymatic cleaning solution with a rust inhibitor, and deionized or distilled water in an ultrasonic cleaner to clean instruments. Use detergents and disinfectants suitable for use with medical and/or dental instruments. Follow manufacturer's instructions related to concentration and length of exposure. Use the processing time recommended by the manufacturer of the detergent and/or ultrasonic system. Ultrasonic solution should be changed daily. Hinged instruments should be cleaned and sterilized in the open position.
- Typical procedure includes: Soak instruments, for the recommended soaking time, in cleaning solution, making sure the instruments are sufficiently immersed. Note: There should not be any contact between the instruments. If applicable, rinse all lumens of the instrument 5 times at the beginning and at the end of the soaking time by application of a single-use syringe (minimum volume 50 ml). If you are using a cassette system, the ultrasonic cleaning time has to be at least 16 minutes, unless a longer exposure time is required by the manufacturer of the detergent. Do not overload the cleaning unit. Remove the instruments from the cleaning solution and post rinse them intensively with deionized or distilled water. Inspect the instruments for a good cleaning result. Dry joint and tip thoroughly with compressed air or towel ensuring joint and tip are free from moisture.

Inspection

- Inspect all instruments after the cleaning and rinsing step for corrosion, damaged surfaces, and impurities. Do not further use damaged instruments. If instruments are still dirty, clean them again.

Sterilization

- CFI Opening Tool is designed for use with steam autoclave.
- Sterilization pouches or cassette systems are recommended.
- Use deionized or distilled water for sterilization.
- The steam sterilization procedure is as follows:

Type of Sterilizer	Temperature	Minimum Time	Drying Time
Gravity Displacement	121°C (250°F)	30 minutes	30 minutes
Dynamic-Air-Removal (Pre-Vacuum)	132°C (270°F)	4 minutes	30 minutes

Storage

- After sterilization, store the instruments in a dry and dust free place, in the clean section of the instrument processing area. Sterilization can only be maintained if the instruments remain packaged or wrapped – impermeable to micro-organisms – following validated standards. The status of the sterilization has to be clearly indicated on the wrapped packages or the containers.

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mdI Europa GmbH Langenhagener Str. 71