## Instructions for Use

### TransForce<sup>®2</sup> Sagittal Arch Developing Appliance

**REF 424-700, 424-728 – 424-740**

**Contre-indications** – This product contains nickel and chromium and should not be used for individuals with known allergic sensitivity to these metals.

TransForce<sup>®2</sup> Sagittal Arch Developing Appliances are manufactured by Ortho Organizers<sup>®</sup> as a series of pre-adjusted fixed appliances designed to correct arch form for patients with contracted dental arches. The appliances deliver smooth continuous force with a pre-activated expansion module that incorporates a nickel titanium coil spring enclosed in a tube.

The TransForce<sup>®2</sup> Sagittal Appliance is designed for antero-posterior arch development in upper or lower arches. The appliance operates on the slide principle and may be used unilaterally or bilaterally to extend arch length. It incorporates bilateral expansion modules, and extends mesially from the molar tube at the gingival level to engage the anterior segment of the dental arch. The expansion module lengthens the arch by using the molars as anchorage while providing force toward the anterior section of the arch. Force of the TransForce<sup>®2</sup> Sagittal Appliance is calibrated to approximately 200 grams of force per expansion module, in accordance with force requirements needed in arch development for Sagittal expansion. No activation is required after fitting the appliance. As a fixed appliance, TransForce<sup>®2</sup> Sagittal Arch Developing Appliances eliminate the challenges with non-compliant patients.


Figure 1—Schematics of TransForce<sup>®2</sup> Sagittal Appliance

#### Arch Development Planner

An Arch Development Planner is provided in each kit to assist in selecting the correct size. This clear template shows how the appliance to scale in both fully compressed and fully expanded forms. A millimeter scale is provided on the template to help you measure.

#### Selecting the Correct Size

1) Measure the existing antero-posterior arch length on each side of the arch by measuring the distance from the molar point, at the mesio-lingual cusp of the molar, to the incisal point, at the midpoint of the gingival papilla between the central incisors.

2) Place planner over the patient's study model to find the size to fit your selected case. Find the "FULLY COMPRESSED" appliance size that best fits inside the lingual outline of the teeth (Figure 2). Compare to the "FULLY EXPANDED" image to see the total amount of expansion capable by that appliance.

Figure 2

For the amount of expansion obtainable with the TransForce<sup>®2</sup> Sagittal Appliance please reference Table 1.

Each appliance is color-coded for easy identification (Table 1).

Catalog No.	Color-Code	Antero-posterior Fully Compressed	Antero-posterior Fully Expanded	Maximum Amount of Expansion
424-728	Green	24 mm	28 mm	4 mm
424-730	Purple	26 mm	30 mm	4 mm
424-732	Pink	28 mm	32 mm	4 mm
424-734	Blue	30 mm	34 mm	4 mm
424-736	Black	30 mm	36 mm	6 mm
424-738	Yellow	32 mm	38 mm	6 mm
424-740	White	34 mm	40 mm	6 mm

Table 1—Measurement for TransForce<sup>®2</sup> Sagittal Appliance

3) The sizes in Table 1 refer to the maximum activation when the appliance is fully compressed compared to fully expanded. In many cases the appliance is not equally compressed on both sides when it is fitted. After activation the appliance equalizes the spaces from molars to incisors on each side. This is an excellent feature for control of dental asymmetry.

4) The larger sizes with 6 mm range of action are normally used in the upper arch by advancing retroclined incisors to accommodate blocked out canines or premolars. This can resolve up to 12 mm of crowding from fully compressed to fully expanded.

5) The smaller sizes with 4 mm range of action are more appropriate in the lower arch to accommodate up to 8 mm crowding by advancing lingually displaced lower incisors. If 4 mm advancement of lower incisors is not indicated it is important to choose the next smaller size to give 2 mm activation to accommodate up to 4 mm of crowding.

6) At any stage in treatment you can de-activate the appliance by crimping the expansion module tubes which will compress them on the wire. If the correct size of appliance is selected, it is not necessary to de-activate the appliance after treatment. It can be used as a passive retainer.

#### Appliance Fitting

1) Within three days prior to placing the appliance, prepare the patient for molar bands by placing separators.

2) On the day of placement, select molar bands with lingual sheaths and confirm correct size of bands.

3) Insert the appliance into both lingual sheath assemblies until the distal extension blade snaps into the latching indent of each sheath (Figure 3). The lingual wires should be occlusal and the expansion module wire should be towards the gingival. Weingart Pliers (Catalog No. 205-202) can be used to insert the placement blade into the sheath.

Figure 3

Inserting the Appliance in the Lingual Sheath

4) For additional support, tie in the appliance by attaching an elastic ligature from the front mesial hook of the blade around to the back of the blade (Figure 4).

Figure 4

Securing the Appliance in the Lingual Sheath Using Elastic Ligature

5) The appliance should be tried on patient's model prior to cementing to verify size and to make any adjustments where necessary (Figure 5).

Figure 5

Fitting the Appliance Prior to Cementing

**Note:** It is easier to attach the molar bands to the TransForce<sup>®2</sup> Sagittal Appliance and fit the appliance in one piece, rather than fitting the bands first then inserting the appliance in the lingual sheaths.

6) Minor adjustments may be required to adapt the appliance to each individual case. The lingual wire can be positioned to help align lingually placed cuspids or bicuspids. You can bend the expansion module wire towards the occlusal margin to prevent tissue impingement.

7) Mix band cement and place in both molar bands.

8) Place and seat one molar band and then the other band. Compress the spring on that side and rotate the anterior wire into position lingual to the anterior teeth. The band on the other side is now displaced distal to the molar. Compress the spring on this side to move the band over the molar and seat the other molar band. Check for secure fit until the cement has hardened. Clean off excess cement after cement has set. The expansion module should be positioned toward the gingival.

9) The device is activated immediately and will provide slow, continuous force to lengthen the arch up to 2 mm every 6 to 8 weeks.

10) After placement and throughout treatment, additional minor adjustments can be made to the appliance intra-orally, if necessary.

#### Maintenance

After fitting the appliance, maintenance should be minimal and may be limited to routine visits of 6- to 8-week intervals to check progress.

#### Retention / Retrait

The appliance can be de-activated and left in place as a retainer after the necessary expansion has been accomplished. To de-activate appliance, crimp the expansion module tube onto the wire. If the correct size of appliance is selected, it is not necessary to de-activate the appliance after treatment. It can be used as a passive retainer.

To remove the appliance, simply stretch and remove elastic ligature, then snap the appliance out of each molar band with Weingart Pliers (Catalog No. 205-202).

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Ortho Organizers and TransForce are registered trademarks of Ortho Organizers, Inc.
Ortho Organizers, Inc.
mdi Europa GmbH