

MASEL®



Doublegard® Mouthguards cushion the shocks of contact sports, guarding against painful labial lacerations.

DOUBLEGARD® THE ORIGINAL PROTECTIVE MOUTHGUARD



When it comes to athletics, mouthguards are more than a matter of looks. They allow coaches and officials to see at a glance whether a player is properly protected. MASEL®'s *Doublegard* Mouthguards conform to NCAA guidelines for colored mouthguards, which is mandatory in most states. The *Doublegard* Mouthguard is also a comfortable, ready-to-wear mouthguard for athletes who wear braces. It cushions the shock of contact sports, guarding against painful labial lacerations.

Doublegard protects your active orthodontic patients


- Stays comfortably in place, ensuring protection
- More comfortable than rubbery "boxer type" mouthguards
- Adapts to teeth and appliances
- Prevents breathing and gagging problems
- Can't be swallowed or dislodged on impact
- Latex-free
- Optional heavy duty strap connects mouthguard to facemask of athletic helmet
- Standard fits ages 12 and older
- Junior fits children ages 11 and under
- Individually wrapped
- 12 per pack

To order by phone (in the U.S.)
Toll-free: 888 851 0533
Fax: 800 888 7244

To order by phone (outside the U.S.)
Direct: + (1) 215 458 9140
Fax: + (1) 215 458 9155

To order online
Email: Sales@Masel.net
Website: Order.MaselOrtho.com

Visit our Website at MaselOrtho.com.



Unlike rubbery mouthguards, Doublegard Mouthguards are soft but durable, comfortable yet tough, and extremely bite resistant

Open channel design adapts easily to teeth and appliances without forcing retrograde tooth movement like impression-type mouthguards may do

Ribbed cushion design

Optional heavy duty strap connects mouthguard to face mask of athletic helmet

Separate upper and lower arch protection lets patients breathe freely

Spring Action Comfort Hinge™ works with the temporomandibular joint—it stays in place, and can't be swallowed or dislodged upon impact

Doublegard with Strap

DESCRIPTION		STANDARD ITEM NO.	JUNIOR ITEM NO.
STANDARD	Translucent	5100-020	5100-040
	Yellow	5100-070	5100-080
	Red	5100-071	5100-081
	Green	5100-072	5100-082
	Royal Blue	5100-073	5100-083
	Black	5100-074	5100-084
	Assortment, 2 of each color	5100-090	5100-091
SUNBURST	Assortment, 3 of each color (Pink, Yellow, Orange, and Green)	5100-100	5100-101
	Assortment, 6 of each color (Gold and Silver)	5100-102	5100-103

Doublegard without Strap

DESCRIPTION		STANDARD ITEM NO.	JUNIOR ITEM NO.
STANDARD	Translucent	5100-010	5100-030
	Yellow	5100-096	—
	Red	5100-092	—
	Green	5100-093	—
	Royal Blue	5100-094	—
	Black	5100-095	—
	Assortment, 2 of each color	5100-097	—

Standard Colors



Sunburst Colors



Metallic Colors

