




Elite[®] Opti-MIM[®]


Mini-Twin[®] Bracket Systems

Imitates a natural tooth shape for superior performance.



Achieve the Ultimate in Control and Comfort

Elite Opti-MIM Mini-Twin Brackets feature a patented tooth-shaped design resulting in a smaller bracket, which delivers uncompromised strength and control throughout the treatment process. Each fully programmed bracket has a compound-contoured base for better bracket to tooth fit. All brackets are manufactured with our *Nickel-Lite*® Stainless Steel (Cobalt Chromium Alloy) material to protect nickel sensitive patients.



Single-Piece Durable Construction

The contour of the bracket allows for an optimum fit on the tooth

Color-coded ID System

For easy identification

Effortless Archwire Changes

Beveled edges of archwire slot provides ease of placement and removal of archwires

Accurate Placement

Distal tie wings are parallel to distal tooth anatomy

Accurate Placement

Mesial tie wings are parallel to mesial anatomy

Precise Vertical Scribe Line


Allows for easy alignment with the long axis of the clinical crown

True Straight Arch System

Advanced engineering is incorporated into each bracket and buccal tube

Elite Opti-MIM Mini-Twin One-piece Bracket offers unsurpassed quality and ease-of-use:

- Available in Roth, Nanda and McLaughlin, Bennett and Trevisi prescriptions*
- Torque-in-base is built into the base of each bracket for an accurate fit of each individual tooth
- Each *Elite Opti-MIM* Bracket and buccal tube feature our patented micro-etched *Grip-LOK*® Base for greater bond strength and reduced bond failure
- Both the regular and low profile *Elite Opti-MIM* Brackets are radiused and polished for maximum patient comfort
- Mesial and distal tie wings are parallel to the tooth's anatomy giving the clinician an additional reference point for accurate placement
- Low profile first and second mini molar tubes feature a stronger bond strength and reduced bond failure**



SEM of *Elite Opti-MIM* Grip-LOK Base




Upper 1st Molar Non-Convertible Buccal Tube

(5x5) All Sets

RX	.018	.022	DESCRIPTION
Roth	702-297	702-397	No hooks - Regular Profile
	702-297LP	702-397LP	No hooks - Low Profile
	702-298	702-398	Hook 3s - Regular Profile
	702-298LP	702-398LP	Hook 3s - Low Profile
	702-299	702-399	Hook 3,4,5s - Regular Profile
	702-299LP	702-399LP	Hook 3,4,5s - Low Profile
McLaughlin, Bennett and Trevisi*	702-299MC	702-399MC	Hook 3,4,5s - Regular Profile
Nanda	702-298RN	702-398RN	Hook 3s - Regular Profile
	702-299RN	702-399RN	Hook 3,4,5s - Regular Profile

Note: Contact customer service for additional prescriptions. All prescription sets feature slim-pad lower bicuspid. *Ortho Organizers*® version of Roth or McLaughlin, Bennett and Trevisi prescriptions is not implied to be an exact version of another system.

TOOTH #	PERMANENT ID FOR REGULAR & LOW PROFILE	COLOR-CODE
U1, U2, U3	Recessed ID on Distal - Gingival Tie Wing	Distal - Gingival Tie Wing
U4HK, U5HK	No Permanent ID	Distal - Gingival Tie Wing
U4, U5	No Permanent ID	Distal - Gingival Tie Wings (Universal)
L1, L2	No Permanent ID	Both - Gingival Tie Wings (Universal)
L3, L4HK(-17T)	Raised ID on Distal - Gingival Tie Wing	Distal - Gingival Tie Wing
L4(17T or -12T)	Raised ID on Distal - Gingival Tie Wing	Both - Gingival Tie Wings (Universal)
L4HK(-12T)	Raised ID on both Gingival Tie Wings	Distal - Gingival Tie Wing
L5HK	Raised ID on both Gingival Tie Wings/Scribe Lines on both Occlusal Tie Wings	Distal - Gingival Tie Wing
L5	Raised ID on both Gingival Tie Wings/Scribe Lines on both Occlusal Tie Wings	Both - Gingival Tie Wings (Universal)












**Ask about our convenient chairside
Quick Pac Bracket Sets**

Elite Opti-MIM

Mini-Twin Brackets









All Prescriptions

Regular Profile

MAXILLARY		COLOR-CODE	TOOTH #	TORQUE	ANGULATION	ITEM NUMBERS			
						.018 R	.018 L	.022 R	.022 L
Centrals		black	U1	+12°	+5°	702-207	702-208	702-307	702-308
				+17°	+4°	702-209	702-210	702-309	702-310
Laterals		pink	U2	+8°	+9°	702-223	702-224	702-323	702-324
				+11°	+9°	702-225	702-226	702-325	702-326
				+10°	+8°	702-227	702-228	702-327	702-328
Cuspids		green	U3	-2°	+11°	702-251	702-252	702-351	702-352
				-2°	+11°	702-251HK	702-252HK	702-351HK	702-352HK
				-7°	+11°	702-249HK	702-250HK	702-349HK	702-350HK
				-7°	+10°	702-253HK	702-254HK	702-353HK	702-354HK
				-7°	+8°	702-255HK	702-256HK	702-355HK	702-356HK
Bicuspid		purple	U4&5s	-7°	0°	702-259	702-260	702-359	702-360
				-7°	0°	702-259HK	702-260HK	702-359HK	702-360HK
MANDIBULAR						.018 R	.018 L	.022 R	.022 L
Anteriors		yellow	L1&2s	0°	0°	702-235	702-236	702-335	702-336
				-6°	0°	702-237	702-238	702-337	702-338
Cuspids		blue	L3	-11°	+5°	702-267	702-268	702-367	702-368
				-11°	+5°	702-267HK	702-268HK	702-367HK	702-368HK
				-6°	+3°	702-271HK	702-272HK	702-371HK	702-372HK
				-6°	+6°	702-273HK	702-274HK	702-373HK	702-374HK
Bicuspid		gray	L4&5s	-12°	0°	702-287SP	702-288SP	702-387SP	702-388SP
				-12°	0°	702-287HKSP	702-288HKSP	702-387HKSP	702-388HKSP
Bicuspid		white	L4&5s	-17°	0°	702-289SP	702-290SP	702-389SP	702-390SP
				-17°	0°	702-289HKSP	702-290HKSP	702-389HKSP	702-390HKSP
Bicuspid		red	L4&5s	-22°	0°	702-293SP	702-294SP	702-393SP	702-394SP
				-22°	0°	702-293HKSP	702-294HKSP	702-393HKSP	702-394HKSP









Note: HK designates bracket comes with a hook.

Selection by Roth Prescription* Roth Regular Profile

MAXILLARY	COLOR-CODE	TOOTH #	TORQUE	ANGULATION	ITEM NUMBERS			
					.018 R	.018 L	.022 R	.022 L
Centrals 	black	U1	+12°	+5°	702-207	702-208	702-307	702-308
Laterals 	pink	U2	+8°	+9°	702-223	702-224	702-323	702-324
Cuspids 	green	U3	-2°	+11°	702-251	702-252	702-351	702-352
		U3	-2°	+11°	702-251HK	702-252HK	702-351HK	702-352HK
Bicuspids 	purple	U4&5s	-7°	0°	702-259	702-260	702-359	702-360
MANDIBULAR					.018 R	.018 L	.022 R	.022 L
Anteriors 	yellow	L1&2s	0°	0°	702-235	702-236	702-335	702-336
Cuspids 	blue	L3	-11°	+5°	702-269	702-270	702-369	702-370
			-11°	+5°	702-267HK	702-268HK	702-367HK	702-368HK
1 st Bicuspids 	white	L4	-17°	0°	702-289SP	702-290SP	702-389SP	702-390SP
			-17°	0°	702-290HKSP	702-289HKSP	702-390HKSP	702-389HKSP
2 nd Bicuspids 	red	L5	-22°	0°	702-294SP	702-293SP	702-394SP	702-393SP
			-22°	0°	702-294HKSP	702-293HKSP	702-394HKSP	702-393HKSP

Note: HK designates bracket comes with a hook.









Selection by Roth Prescription* Roth Low Profile

MAXILLARY	COLOR-CODE	TOOTH #	TORQUE	ANGULATION	ITEM NUMBERS			
					.018 R	.018 L	.022 R	.022 L
Centrals 	black	U1	+12°	+5°	702-205	702-206	702-305	702-306
Laterals 	pink	U2	+8°	+9°	702-221	702-222	702-321	702-322
Cuspids 	green	U3	-2°	+11°	702-251	702-252	702-351	702-352
			-2°	+11°	702-251HK	702-252HK	702-351HK	702-352HK
Bicuspids 	purple	U4&5s	-7°	0°	702-259	702-260	702-359	702-360
		U4&5s	-7°	0°	702-259HK	702-260HK	702-359HK	702-360HK
MANDIBULAR					.018 R	.018 L	.022 R	.022 L
Anteriors 	yellow	L1&2s	0°	0°	702-233	702-234	702-333	702-334
Cuspids 	blue	L3	-11°	+5°	702-269	702-270	702-369	702-370
			-11°	+5°	702-269HK	702-270HK	702-369HK	702-370HK
1 st Bicuspids 	white	L4	-17°	0°	702-289SP	702-290SP	702-389SP	702-390SP
			-17°	0°	702-289HKSP	702-290HKSP	702-389HKSP	702-390HKSP
2 nd Bicuspids 	red	L5	-22°	0°	702-293SP	702-294SP	702-393SP	702-394SP
			-22°	0°	702-293HKSP	702-294HKSP	702-393HKSP	702-394HKSP

Note: Centrals and Laterals are available in a lower profile and can be used with any prescription of our *Elite Opti-MIM* regular profile cuspids and bicuspid brackets. For cuspids and bicuspids refer to Regular Profile Roth Prescription. HK designates bracket comes with a hook.

Selection by McLaughlin, Bennett and Trevisi Prescription*









MBT Regular Profile

MAXILLARY	COLOR-CODE	TOOTH #	TORQUE	ANGULATION	ITEM NUMBERS				
					.018 R	.018 L	.022 R	.022 L	
Centrals		black	U1	+17°	+4°	702-209	702-210	702-309	702-310
Laterals		pink	U2	+10°	+8°	702-227	702-228	702-327	702-328
Cuspids		green	U3	-7°	+8°	702-255HK	702-256HK	702-355HK	702-356HK
Bicuspids		purple	U4&5s	-7°	0°	702-259 702-259HK	702-260 702-260HK	702-359 702-359HK	702-360 702-360HK
MANDIBULAR						.018 R	.018 L	.022 R	.022 L
Anteriors		yellow	L1&2s	-6°	0°	702-237	702-238	702-337	702-338
Cuspids		blue	L3	-6°	+3°	702-271HK	702-272HK	702-371HK	702-372HK
1 st Bicuspids		gray	L4	-12°	0°	702-287SP 702-287HKSP	702-288SP 702-288HKSP	702-387SP 702-387HKSP	702-388SP 702-388HKSP
2 nd Bicuspids		red	L5	-17°	0°	702-289SP 702-289HKSP	702-290SP 702-290HKSP	702-389SP 702-389HKSP	702-390SP 702-390HKSP

Note: HK designates bracket comes with a hook.

Selection by Nanda Prescription*





Nanda Regular Profile

MAXILLARY	COLOR-CODE	TOOTH #	TORQUE	ANGULATION	ITEM NUMBERS				
					.018 R	.018 L	.022 R	.022 L	
Centrals		black	U1	+17°	+4°	702-209	702-210	702-309	702-310
Laterals		pink	U2	+11°	+9°	702-225	702-226	702-325	702-326
Cuspids		green	U3	-7°	+10°	702-253HK	702-254HK	702-353HK	702-354HK
Bicuspids		purple	U4&5s	-7°	0°	702-259	702-260	702-359	702-360
			U4&5s	-7°	0°	702-259HK	702-260HK	702-359HK	702-360HK
MANDIBULAR						.018 R	.018 L	.022 R	.022 L
Anteriors		yellow	L1&2s	0°	0°	702-235	702-236	702-335	702-336
Cuspids		blue	L3	-11°	+5°	702-267HK	702-268HK	702-367HK	702-368HK
1 st Bicuspids		white	L4	-17°	0°	702-289SP	702-290SP	702-389SP	702-390SP
			L4	-17°	0°	702-289HKSP	702-290HKSP	702-389HKSP	702-390HKSP
2 nd Bicuspids		red	L5	-22°	0°	702-293SP	702-294SP	702-393SP	702-394SP
			L5	-22°	0°	702-293HKSP	702-294HKSP	702-393HKSP	702-394HKSP

Note: HK designates bracket comes with a hook.



Elite Opti-MIM Direct Bond Molar Tubes – Non Convertible

All *Elite Opti-MIM* Molar Tubes are made from *Nickel-Lite* Material and have a patented *Grip-LOK* Base, improving bond strength and greatly reducing bond failure. The defined buccal indent assists with accurate placement and stability.

MAXILLARY	TOOTH #	TORQUE	ANGULATION	OFFSET	ITEM NUMBERS				
					.018 R	.018 L	.022 R	.022 L	
1 st Molar		U6	-10°	0°	+10°	901-401	901-402	901-451	901-562
		U6	-10°	0°	+10°	901-403	901-404	901-453	901-454
2 nd Molar		U7	-14°	0°	+14°	951-205	951-206	951-255	951-256
		U7	-14°	0°	+14°	951-207	951-208	951-257	951-258
Edgewise		Tip U7	Torque 0°	0°	+10°	951-201	951-202	951-251	951-252
MANDIBULAR						.018 R	.018 L	.022 R	.022 L
1 st Molar		L6	-12°	0°	+4°	901-411	901-412	901-461	901-462
		L6	-25°	0°	+4°	901-413	901-414	901-463	901-464
2 nd Molar		L7	-12°	0°	+4°	951-209	951-210	951-259	951-260
		L7	-25°	0°	+4°	951-211	951-212	951-261	951-262
Edgewise		Tip L7	Torque 0°	0°	+4°	951-203	951-204	951-253	951-254



Elite Opti-MIM Direct Bond Micro Molar Tubes

Ortho Organizers' Direct Bond Micro Molar Tubes have a reduced pad area enhanced with our micro-etched *Grip-LOK* Base. The defined buccal indent assists with accurate placement and stability. The micro tubes may be used on first and second molars. The shorter flanges allow the buccal tube a more precise interface with the tooth.

MAXILLARY	TOOTH #	TORQUE	ANGULATION	OFFSET	ITEM NUMBERS				
					.018 R	.018 L	.022 R	.022 L	
1 st Molar/ 2 nd Molar		U6/7	-10°	0°	+10°	901-401MP	901-402MP	901-451MP	901-452MP
		U6/7	-14°	0°	+14°	901-403MP	901-404MP	901-453MP	901-454MP
MANDIBULAR						.018 R	.018 L	.022 R	.022 L
1 st Molar/ 2 nd Molar		L6/7	-12°	0°	+4°	901-411MP	901-412MP	901-461MP	901-462MP
		L6/7	-25°	0°	+4°	901-413MP	901-414MP	901-463MP	901-464MP

Elite Opti-MIM Direct Bond Solo Tubes

One-piece, one cusp solo tubes are about half the length of a traditional buccal tube, and are designed to be bonded to the mesial cusp of the 1st or 2nd molar. The compound contour of the solo tube and the patented micro-etched *Grip-LOK* Base creates a stable fit for placement.

MAXILLARY	TOOTH #	TORQUE	ANGULATION	OFFSET	ITEM NUMBERS				
					.018 R	.018 L	.022 R	.022 L	
1 st Molar/ 2 nd Molar		U6/7	0°	0°	0°	901-501	901-502	901-551	901-552
		U6/7	-10°	0°	0°	901-503	901-504	901-553	901-554
MANDIBULAR						.018 R	.018 L	.022 R	.022 L
1 st Molar/ 2 nd Molar		L6/7	0°	0°	0°	901-511	901-512	901-561	901-562
		L6/7	-12°	0°	0°	901-513	901-514	901-563	901-564
			-25°	0°	0°	901-515	901-516	901-565	901-566

Please contact *Ortho Organizers* for a complete listing of our buccal tube product selections which includes convertible or non-convertible, double or triple, standard, lip bumper or facebow. For added convenience, all selections are available loose, direct bond or pre-welded to our bands.



Subscribe to our eNews: OrthoOrganizers.com/eNews

For more information on our products, orthodontic laboratory services, and educational offerings, please contact us:

Customer Service: 888.851.0533

In the U.S.: 800.547.2000

Outside the U.S.: +(1) 760 448 8600

Canada: CERUM 800.661.9567

To fax an order: 800.888.7244

To fax an order outside the U.S.: +(1) 760 448 8607

To email an order: USASales@OrthoOrganizers.com

To email an order outside the U.S.: IntlOrders@OrthoOrganizers.com



ORTHO[®]
Organizers

A  HENRY SCHEIN[®] COMPANY

1822 Aston Avenue
Carlsbad, CA 92008
USA
OrthoOrganizers.com

* Ortho Organizers prescription does not imply endorsement by the doctor(s).
** Study available upon request.

