

# ANDREWS<sup>2</sup>™ ARCHWIRES

## ORTHODONTIC SCIENCE REALIZED.

Optimal arches are essential for optimal occlusion. *Andrews<sup>2</sup> Archwires* have the features needed to treat arches to optimal standards efficiently.

- .014, .016 *Nitanium*® Super Elastic Wires have flexibility for initial aligning and rotations
- .014, .016 stainless steel wires have a balance of flexibility and rigidity for initial aligning and leveling
- .016, .018 stainless steel wires have a balance of flexibility and rigidity for leveling, correcting arch shape, and for mesiodistal translation
- .018 x .025 Ortholoy-Blue (heat treatable) stainless steel wires are for inclining and fine tuning

### Flexibility:

Arches can be effectively and efficiently treated to optimal standards by using the most appropriate *Andrews<sup>2</sup> Appliance* prescription in combination with the most suitable *Andrews<sup>2</sup> Archwires* in terms of flexibility.

### Shape:




The *Andrews<sup>2</sup> Archwires* are preformed to match the average shape of the 120 non-orthodontic normal study casts used in the "Six Keys to Normal (Optimal) Occlusion" study. Like any preformed archwire, its shape must be customized to each patient's unique anatomy.

### Length:

The overall length can usually be customized with one cut.

### Pre-Posted Archwire Option:

For each arch, the small, medium, and large distances between the posts means one of the 3 will fit 90% of the normal patient population.\*





# ANDREWS<sup>2</sup> ARCHWIRES

Archwire	NITANIUM SUPER ELASTIC ARCHWIRE	STAINLESS STEEL ARCHWIRE	ORTHOLOY-BLUE ARCHWIRE
Arch Form	Andrews <sup>2</sup> Arch Form	Andrews <sup>2</sup> Arch Form	Andrews <sup>2</sup> Arch Form
Packaged	Single Pack 10/Pack	Single Pack 10/Pack	Single Pack 10/Pack
.014 Upper	112-140	112-100	—
.014 Lower	112-141	112-101	—
.016 Upper	112-142	112-102	—
.016 Lower	112-143	112-103	—
.018 Upper	112-144	112-104	—
.018 Lower	112-145	112-105	—
.017 x .025 Upper**	112-146	—	—
.017 x .025 Lower**	112-147	—	—
.018 x .025 Upper	—	—	112-172
.018 x .025 Lower	—	—	112-173

Archwire	ORTHOLOY-BLUE PRE-POSTED ARCHWIRE
Arch Form	Andrews <sup>2</sup> Arch Form
Packaged	Single Pack 10/Pack
.018 x .025 Pre-Posted Upper	
33 mm posted	112-174
37 mm posted	112-178
41 mm posted	112-182
.018 x .025 Pre-Posted Lower	
24 mm posted	112-175
28 mm posted	112-179
32 mm posted	112-183

Archwire	STAINLESS STEEL PRE-POSTED ARCHWIRE
Arch Form	Andrews <sup>2</sup> Arch Form
Packaged	Single Pack 10/Pack
.016 Pre-Posted Upper	
33 mm posted	112-106
37 mm posted	112-116
41 mm posted	112-126
.016 Pre-Posted Lower	
24 mm posted	112-107
28 mm posted	112-117
32 mm posted	112-127
.018 Pre-Posted Upper	
33 mm posted	112-108
37 mm posted	112-118
41 mm posted	112-128
.018 Pre-Posted Lower	
24 mm posted	112-109
28 mm posted	112-119
32 mm posted	112-129



Andrews<sup>2</sup> Archwires are great for use with the Andrews<sup>2</sup> Appliance.

**Andrews<sup>2</sup>  
Appliance**  
orthodontic science realized.

For more information on Ortho Organizers' products and education offerings, please contact us:  
In the U.S. 888.851.0533 | Outside the U.S. +(1) 760 448 8600 | Canada: CERUM 800.661.9567  
To fax an order: 800.888.7244 | To email an order: [USASales@OrthoOrganizers.com](mailto:USASales@OrthoOrganizers.com)



Subscribe to our eNews: [OrthoOrganizers.com/eNews](http://OrthoOrganizers.com/eNews)

**ORTHOORGANIZERS.COM**

\*Call for availability.

\*\* To be used for short periods for aligning slots.

© 2013 Ortho Organizers, Inc. All rights reserved. PN 999-129 05/13

**ORTHO<sup>2</sup>**  
Organizers  
A HENRY SCHEIN® COMPANY

