



Nitanium Molar Rotator²TM

Multi-Purpose Finishing Appliance

For patients requiring upper first molar rotation and stabilization

The *Nitanium Molar Rotator*² Appliance provides predictable rotation, torquing and expansion in a pre-programmed appliance.

Advances in thermally activated nickel titanium provide the activating forces in this truly unique appliance. After insertion into the maxillary first molar lingual sheaths, the appliance provides a constant force to ensure the desirable end result.

The *Nitanium Molar Rotator*² Appliance offers many benefits, including the following:

- Slower expansion has shown to lower incidents of relapse*
- Thermally-activated *Nitanium*[®] for total control eliminating patient's compliance
- Ease of placement and removal by the doctor
- Orthology arm and loops provide adjustability for fine-tuning and finishing
- No laboratory work is required – available in 10 sizes
- Excellent for mixed dentition and adolescent cases
- Consistent and predictable end results

Nitanium Molar Rotator² Appliance

Individual Rotators *, ** (1 per pack)

26 mm Rotator	•	101-781
28 mm Rotator	•	101-782
30 mm Rotator	•	101-783
32 mm Rotator	•	101-784
34 mm Rotator	•	101-785
36 mm Rotator	•	101-786
38 mm Rotator	•	101-787
40 mm Rotator	•	101-788
42 mm Rotator	•	101-789
44 mm Rotator	•	101-790

* For easy identification, each rotator has size numbers permanently stamped on the assembly.

** Rotators are for single-use only.

Nitanium Molar Rotator² Appliance Full Range Kit • 101-780

1 Rotator of each of the following 10 sizes:
26 mm, 28 mm, 30 mm, 32 mm, 34 mm,
36 mm, 38 mm, 40 mm, 42 mm, and 44 mm
1 can of Ortho Freez (10 oz)
1 acrylic storage holder
Instruction manual

Individual Components

Ortho Freez (10 oz can) • 100-788
Rotator Adjustment Plier • 200-450

Auxiliary Products

The **Nitanium Molar Rotator²** Appliance simply slips into any of the following weldable Lingual Sheaths.



971-032

.036 slotted, horizontal,
latching indent.
Fits left or right.
(10 per pack)



UR/LL • 971-034 UL/LR • 971-035

.036 slotted, horizontal,
latching indent hook.
(10 per pack)



UR/LL • 971-036 UL/LR • 971-037

.036 horizontal hook,
slotted, latching indent.
(10 per pack)



UR/LL • 971-038 UL/LR • 971-039

.036 vertical hooks.
(10 per pack)

Individual band kits are available with a welded Lingual Sheath and Buccal Tube.



Subscribe to our eNews: OrthoOrganizers.com/eNews

For more information on our products, orthodontic laboratory services, and educational offerings, please contact us:

Customer Service: 888.851.0533

In the U.S.: 800.547.2000

Outside the U.S.: +(1) 760 448 8600

Canada: CERUM 800.661.9567

To fax an order: 800.888.7244

To fax an order outside the U.S.: +(1) 760 448 8607

To email an order: USASales@OrthoOrganizers.com

To email an order outside the U.S.: IntlOrders@OrthoOrganizers.com



ORTHO
Organizers

A HENRY SCHEIN® COMPANY

1822 Aston Avenue
Carlsbad, CA 92008
USA
OrthoOrganizers.com

