



ADHESIVES/CEMENTS | F







| | |
|------------------------|------|
| Adhesives | F-2 |
| Band Cements | F-7 |
| Curing Light | F-9 |
| Auxiliaries | F-9 |
| Conditioners | F-11 |

| ITEM | DESCRIPTION | CATALOG NO. |
|---|--|-------------|
|  | <p>Everest™ Light-Cure Orthodontic Adhesive</p> <p>The smooth, tacky composite resin delivers ideal handling characteristics which eliminates bracket float and allows for easy flash clean up. Adhesive provides the optimum properties for bonding metal, ceramic, or plastic brackets using light-cure technology. Light-activation provides accurate bracket placement ensured by allowing for sufficient working time. Produces a strong light-cured chemical bond which also debonds predictably. Archwires can be placed immediately after fully curing the brackets. Contains fluoride to help protect against potential decalcification.</p> | |
| | <p>Everest Light-Cure Kit in Syringes</p> <p>Kit Includes:</p> <ul style="list-style-type: none"> • 4 Push Syringes (3.5 grams ea.) • 1 Gel Etch Syringe (4 grams) • 1 Sealant Resin (7 ml) • 5 Syringe Tips (22-Gauge) • 2 Mixing Wells (10-Cavity) • 1 Mixing Pad • 50 Brushes • 1 Instruction For Use Sheet | 490-300 |
| | <p>Everest Light-Cure Kit in Capsules</p> <p>Kit Includes:</p> <ul style="list-style-type: none"> • 40 Capsules (.35 grams each) • 1 Sealant Resin (7 ml) • 1 Green Liquid Etchant (9 grams) • 1 Mixing Pad • 50 Brushes • 1 Instruction For Use Sheet | 490-320 |
| | <p>Everest Dispensing Gun</p> <p>1 Capsule Gun</p> | 490-220 |
| | <p>Everest Light-Cure Mini Kit</p> <p>Kit Includes:</p> <ul style="list-style-type: none"> • 1 Push Syringe (3.5 grams) • 1 Sealant Resin (3 ml) • 1 Green Gel Etchant Syringe (4 grams) • 5 Syringe Tips (22-Gauge) • 5 Bonding Brushes • 1 Instruction For Use Sheet | 490-310 |
| | <p>Everest Light-Cure Reorder Components</p> | |
| | <p>Light-Cure Adhesive Syringe Pack (4 Syringes, 3.5 grams ea.)</p> | 490-111 |
| | <p>4 Light-Cure Adhesive Single Syringe (3.5 grams ea.)</p> | 490-112 |
| | <p>40 Light-Cure Single Use Capsules (.35 grams ea.)</p> | 490-221 |
| | <p>1 Sealant Resin (7 ml)</p> | 490-200 |
| <p>1 Gel Etchant Large Bottle (330 grams)</p> | 490-209 | |
| <p>1 Gel Etchant Small Bottle (57 grams)</p> | 490-210 | |
| <p>5 Mixing Pads</p> | 490-505 | |
| <p>1 Brush Holder</p> | 490-504 | |
| <p>100 Disposable Brushes</p> | 490-510 | |
| <p>50 Disposable Long-Handled Brushes</p> | 490-520 | |


| ITEM | DESCRIPTION | CATALOG NO. |
|---|--|--|
|  | <p>Everest No-Mix Orthodontic Adhesive (Primer Activated) Bonding System Provides the optimum properties for bonding ceramic, metal, or plastic brackets using a primer activated no-mix bonding system. This is a chemically-cured system which does not require pre-mixing. The thick, tacky adhesive resin delivers ideal handling characteristics which eliminate bracket flotation and allow for easy flash clean up. Archwires can be placed 5 minutes after bonding the last bracket. Contains fluoride to help protect against potential decalcification.</p> <p>Everest No-Mix Orthodontic Adhesive Kit in Syringes (Primer-Activated) Kit Includes:</p> <ul style="list-style-type: none"> • 4 Push Syringes (3.5 grams ea.) • 1 Enamel Catalyst Primer (9 ml) • 1 Gel Etchant Syringe (4 grams) • 5 Syringe Tips (22-Gauge) • 2 Mixing Wells (10-Cavity) • 1 Mixing Pad • 50 Brushes • 1 Instruction For Use Sheet <hr/> <p>Everest No-Mix Orthodontic Adhesive Mini Kit in Syringe (Primer-Activated) Kit Includes:</p> <ul style="list-style-type: none"> • 1 Push Syringe (3.5 grams) • 1 Enamel Catalyst Primer (3 ml) • 1 Gel Etchant Syringe (4 grams) • 5 Syringe Tips (22-Gauge) • 5 Bonding Brushes • 1 Instruction For Use Sheet | <p>490-340</p> <p>490-350</p> |
| | <p>Everest Etchant Bonding Supplies Gel stays where placed without overflowing and provides highly visible contrast against tooth surfaces. Non-drying water soluble consistency with ideal medium viscosity. Traditional 37% Phosphoric acid.</p> <p>Everest Gel Etch Kit, Economy Size Kit Includes:</p> <ul style="list-style-type: none"> • 4 Syringes (4 grams ea.) • 20 Plastic Flowthru Tips (22-Gauge) <hr/> <p>Everest Gel Etch Mini Kit Kit Includes:</p> <ul style="list-style-type: none"> • 1 Syringe (4 grams) • 5 Plastic Flowthru Tips (22-Gauge) <hr/> <p>Everest Etchant Reorder Components</p> <ul style="list-style-type: none"> 1 Gel Etchant Large Bottle (330 grams) 1 Gel Etchant Small Bottle (57 grams) 1 Liquid Gel Etchant Large Bottle (29.6 ml) 1 Liquid Gel Etchant Small Bottle (7 ml) 25 Etchant Gel Needle Tips | <p>490-360</p> <p>490-370</p> <p>490-209 490-210 490-211 490-212 490-235</p> |


ADHESIVES

| ITEM | DESCRIPTION | CATALOG NO. |
|---|--|--|
|  | <p>Fuji ORTHO™ Automix Light-Cure (LC) Bracket Adhesive with Fluoride The Automix System simplifies application – ensuring a perfect mix each time. With the Automix System there is no need for hand mixing or an amalgamator. The adhesive mixes while it dispenses making application and preparation easy and stress free.</p> <p>Longer working time</p> <ul style="list-style-type: none"> • Five full minutes of working time <p>Always the perfect mix</p> <ul style="list-style-type: none"> • Automix system makes preparation stress free, no measuring or hand mixing needed <p>No white spots</p> <ul style="list-style-type: none"> • Resin-reinforced orthodontic glass-ionomer adhesive - provides a long-term fluoride release which aides in the prevention of white spot lesions <p>Moisture tolerant</p> <ul style="list-style-type: none"> • Wet field compatibility <p>Blue color for bands</p> <ul style="list-style-type: none"> • Makes cement easier to distinguish from the tooth surface for easy removal <hr/> <p>Fuji ORTHO Automix Light-Cure (LC) Starter Kit (Starts approximately 13 cases)</p> <p>Kit Includes:</p> <ul style="list-style-type: none"> • 1 Paste Pak Cartridge (13.3 grams) • 1 Paste Pak Dispenser • 2 Ortho Gel Conditioner Syringes (1.4 grams each) • 20 Automix Mixing Tips • 5 Brushes for Conditioner <hr/> <p>Fuji ORTHO Automix Light-Cure (LC) Set (Starts approximately 13 cases)</p> <p>Kit Includes:</p> <ul style="list-style-type: none"> • 1 Paste Pak Cartridge (13.3 grams) • 2 Ortho Gel Conditioner Syringes (1.4 grams each) • 20 Automix Mixing Tips • 15 Brushes for Conditioner <p>Automix (LC) Reorder Components</p> <p>2 Paste Pak Cartridges (13.3 grams each) 20 Automix Mixing Tips 3 Ortho Gel Conditioner Syringes (1.4 grams each) with 25 Brushes for Conditioner</p> | <p>480-112</p> <p>480-114</p> <p>480-113 480-115 480-116</p> |

| ITEM | DESCRIPTION | CATALOG NO. |
|--|--|---|
|  | <p>Fuji ORTHO Light-Cure (LC) Bonding Adhesives with Fluoride This light-cured resin modified glass-ionomer adhesive is suited for all bonding procedures. Placement in a moist field simplifies application and bonding. Light cures in as little as 40 seconds. Releases fluoride to help prevent decalcification. Pre-capsulation produces ideal bonding consistency in every mix. Bonding procedures can be performed more easily and efficiently.</p> <p>Intro LC Kit (starts approx. 14 cases) Kit Includes:</p> <ul style="list-style-type: none"> • Powder (15 grams) • Liquid (8 grams) • 1 Mixing Pad • 1 Measuring Spoon <hr/> <p>Standard LC Kit (starts approx. 34 cases) Kit Includes:</p> <ul style="list-style-type: none"> • Powder (40 grams) • Liquid (16 grams) • 1 Mixing Pad • 1 Measuring Spoon <hr/> <p>Intro and Standard LC Reorder Components Conditioner (25 grams) Powder (40 grams) Liquid (8 grams)</p> | <p>480-101</p> <p>480-102</p> <p>480-108 480-106 480-107</p> |
|   | <p>LC Capsules Starter Kit (starts approx. 25 cases) Pre-capsulation produces ideal bonding consistency in every mix. Kit Includes:</p> <ul style="list-style-type: none"> • 50 Pre-measured and Sealed Capsules (that can be activated in an amalgamator) • 1 Metal Capsule Applicator • 1 Bottle of Conditioner <p>LC Capsule Reorder Components Capsules Refill Packages (50 capsules/starts approx. 25 cases) Metal Capsule Applicator</p> <p>Fuji ORTHO Self-Cure (SC) Bonding Adhesives with Fluoride (Bond in Moist Field) This self-cured resin modified glass ionomer adhesive is ideally suited for bonding brackets. Placement in a moist field simplifies application and bonding procedures. Fluoride release helps prevent against decalcification and decay. Can be used with an etch or non-etch technique. As a self-curing adhesive, it is particularly useful in bonding situations in which light curing is neither possible nor desired.</p> <p>Fuji ORTHO Standard SC Kit (starts approx. 34 cases) Kit Includes:</p> <ul style="list-style-type: none"> • 1 Bottle of Powder (40 grams) • 2 Bottles of Liquid (8 grams, 6.8 ml each) • 1 Scoop • 1 Mixing Pad <p>Fuji ORTHO SC Reorder Components Powder (40 grams) Liquid (8 grams) Conditioner (25 grams)</p> | <p>480-133</p> <p>480-131</p> <p>480-121</p> <p>480-126 480-127 480-108</p> |








BAND CEMENTS

| ITEM | DESCRIPTION | CATALOG NO. |
|--|---|---|
|  | <p>Everest Light-Cure Band Cement Provides a longer working time and a rapid setting time with its light-activated band cement in a single paste system. Can be placed directly into bands. Excess flash is easy to clean up prior to the light curing and the cement is easy to see due to the blue color. Contains fluoride to help protect against potential decalcification.</p> <p>Everest Light-Cure Orthodontic Band Cement Kit Kit Includes: • 6 Push Syringes (5 grams ea.)</p> <hr/> <p>Everest Light-Cure Orthodontic Band Cement Reorder Components 1 Light-Cure Band Cement Refill Syringe (5 grams)</p> | <p>490-390</p> <p>490-230</p> |
| | <p>Everest Powder/Liquid Band Cement Everest chemically bonds to enamel without the use of etchants. Using a two-part glass ionomer band cement system of powder and liquid, it has a low film thickness and will not wash out under bands. The final setting time is 6 to 7 minutes. Contains fluoride to help protect against potential decalcification.</p> <p>Everest Glass Ionomer Cement Economy Kit Kit Includes: • 1 Powder (100 grams) • 2 Liquid (50 grams ea.) • 1 Mixing Pad • 1 Measuring Scoop • 1 Instruction For Use Sheet</p> <hr/> <p>Everest Glass Ionomer Cement Starter Kit Kit Includes: • 1 Powder (50 grams) • 1 Liquid (50 grams) • 1 Mixing Pad • 1 Measuring Scoop • 1 Instruction For Use Sheet</p> <hr/> <p>Everest Glass Ionomer Cement Reorder Components 1 Glass Ionomer Powder (50 grams) 1 Glass Ionomer Liquid (50 ml)</p> | <p>490-380</p> <p>490-385</p> <p>490-270 490-271</p> |
| | <p>Everest Miscellaneous Reorder Components 50 Disposable Spatulas 5 Brush Holders 5 Mixing Pads 5 Measuring Scoops 100 Disposable Brushes 50 Disposable Long-Handled Brushes 10 Mixing Wells (2-Cavity) 1 Dispensing Gun 25 Etchant Gel Needle Tips & LV Flowable Tips</p> | <p>490-500 490-504 490-505 490-506 490-510 490-520 490-535 490-220 490-235</p> |







| ITEM | DESCRIPTION | CATALOG NO. |
|---|---|--|
|  | <p>Fuji ORTHO Automix Light-Cure Band Cement with Fluoride Makes bonding faster and easier, the Automix systems dispenses the exact ratio so there is a perfect mix every time. Application is stress free with no need to measure, mix, or etch and paste placement is in a wet field.</p> <p>Longer working time</p> <ul style="list-style-type: none"> • Five full minutes of working time <p>Always the perfect mix</p> <ul style="list-style-type: none"> • Automix system makes preparation stress free, no measuring or hand mixing needed <p>No white spots</p> <ul style="list-style-type: none"> • Resin-reinforced orthodontic glass-ionomer adhesive - provides a long-term fluoride release which aides in the prevention of white spot lesions <p>Moisture tolerant</p> <ul style="list-style-type: none"> • Wet field compatibility <p>Blue color for bands</p> <ul style="list-style-type: none"> • Makes cement easier to distinguish from the tooth surface for easy removal <p>Fuji ORTHO Automix Light-Cure Band Cement Starter Kit (enough to cement 56 bands) Kit Includes:</p> <ul style="list-style-type: none"> • 2 Paste Pak Cartridges (13.3 grams each) • 1 Paste Pak Dispenser • 20 Automix Mixing Tips <p>Fuji Ortho Automix Light-Cure Band Cement Reorder Components</p> <ul style="list-style-type: none"> • 2 Paste Pak Cartridges (13.3 grams each) • 1 Mixing Pad | <p>480-110</p> <p>480-111</p> |
|  | <p>Fuji ORTHO Band Cement Self-Cure (SC) Makes bonding faster and easier – the dispenser produces two pastes to the exact ratio so a strong bond can be achieved every time. The doctor or assistant mixes the paste for 15 seconds to a smooth, creamy consistency. The paste stays in place without running and is good for up to 3½ minutes. No etching is required and paste is placed in a wet field. Each kit contains enough paste to cement approximately 62 bands.</p> <p>Fuji ORTHO Band Cement Self-Cure Kit (enough to cement 62 bands) Kit Includes:</p> <ul style="list-style-type: none"> • 2 Paste Pak Cartridges (13.3 grams each) • 1 Paste Pak Dispenser • 1 Mixing Pad • 1 Mixing Spatula <p>Fuji ORTHO Band Cement Self-Cure Reorder Components</p> <p>2 Paste Pak Cartridges (13.3 grams each) 1 Paste Pak Dispenser</p> | <p>480-140</p> <p>480-141 480-142</p> |

CURING LIGHT/AUXILIARIES

| ITEM | DESCRIPTION | CATALOG NO. |
|---|--|--|
|  | <p>Fusion™ 4 Ortho LED Curing Light*</p> <p><i>Fusion 4 Ortho LED Curing Light is a high power, lower dispersion curing light that enables rapid bracket cures in only 3-5 seconds. This modular light is conveniently sized to fit in the palm of your hand and features a rechargeable smart lithium-ion battery for cordless use. This 2700-mW/cm² laser-sharp focused beam is activated through a sleek metal keypad that includes raised buttons including a large ON/OFF button for easy activation.</i></p> <p><i>Fusion Curing Light Kit (Silver)</i></p> <p><i>Fusion Curing Light Kit (Pink)</i></p> <p><i>Fusion Curing Caps-13 mm (50/Pack)</i></p> <p><i>Fusion Barrier Sleeves (100/pack)</i></p> <p><i>Fusion Eye Shield</i></p> <div style="display: flex; justify-content: space-around; align-items: center;">      </div> <p><small>* Features a universal adapter for international use</small></p> | <p>506-0360</p> <p>506-0361</p> <p>506-0375</p> <p>506-0370</p> <p>506-0372</p> |
|  | <p>Mini Sponge Applicators (200 Sponges in Plastic Jar)</p> <p>Specially designed for application of liquids, etches, and resins to teeth. Using tweezers, each sponge can be picked out from premoistened jar and dipped in liquid before applying to desired area. Mini sponges may be used during bonding to apply etchant or any liquid, such as primers, conditioners, etc.</p> | <p>480-015</p> |
|  | <p>Gel Etch Econo-Syringe*</p> <p>Gel etch is made up of a green, 37% phosphoric acid formula and is dispensed from a large 12 gram syringe.</p> <p> * Ground shipment recommended – This product is classified as hazardous material as per the U.S. Department of Transportation Pipeline and Hazardous Materials Safety Administration. Extra shipping charges may apply.</p> | <p>470-162</p> |

| ITEM | DESCRIPTION | CATALOG NO. |
|---|--|--|
|  | <p>Cheek Retractors</p> <p>For treatment with teeth in open or closed bite. Self-retaining. Can be cold sterilized. Excellent for photography.</p> <ul style="list-style-type: none"> • Adult Cheek Retractors (2 per pack) • Teen Cheek Retractors (2 per pack) • Pedo Cheek Retractors (2 per pack) • Assortment Pack (one of each size) | <p>300-020 300-022 300-021 300-023</p> |
|  | <p>Adjustable Cheek Retractors</p> <p>Unique bow design allows adjustment to ensure comfortable fit. Can withstand temperatures up to 400°F. Fully sterilizable. (1 per pack)</p> | <p>300-015</p> |
|  | <p>Lip Retractors</p> <p>Transparent plastic lip retractors have long handles so they can be held easily. Excellent for photography. Rounded edges. Can be cold sterilized. (2 per pack)</p> <ul style="list-style-type: none"> • Standard Lip Retractors • Mini Lip Retractors | <p>300-025 300-026</p> |
|  | <p>Double End Lip Retractors</p> <p>Features both a pedo and teen/adult size. Comfortable and easy to hold, yet allows an excellent field of vision for treatment and photographs. Can withstand temperatures of 370°F/188°C. Amber tint identifies it as an autoclavable plastic instrument. (2 per pack)</p> | <p>300-027</p> |
|  | <p>Retract EEZ® Photo Lip Retraction</p> <p>Excellent for occlusal photographs with easy to grip handles, a great anatomical fit, plus angled handles. Can be cold sterilized, as well as autoclaved, though heat sterilization may cause hazing of the material.</p> | <p>300-029</p> |
|  | <p>Combo Cheek Expander</p> <p>Perform two tasks at once; expand the cheeks while simultaneously depressing the tongue for unrestricted access to the whole mouth. Can withstand temperatures of 370°F/188°C. Amber tint identifies it as an autoclavable plastic instrument.</p> <ul style="list-style-type: none"> • Small • Medium • Large | <p>300-016 300-017 300-018</p> |
|  | <p>Tongue Shield – Latex Free</p> <p>Combination of two bite blocks and a vertical shield. The shield prevents the tongue from touching the lingual surface of the anterior teeth. Center hole allows for insertion of saliva ejector. Excellent during bonding procedures. Can be cold sterilized. (Pink, 3 per pack)</p> <ul style="list-style-type: none"> • Adult/Teen • Child/Pedo | <p>300-030 300-031</p> |

AUXILIARIES/CONDITIONERS

| ITEM | DESCRIPTION | CATALOG NO. |
|---|--|--|
|  | <p>Disposable Tongue Block</p> <p>Available in two sizes to accommodate pedo or teen/adults. These strawberry scented tongue blocks keep the mouth open, with the jaw muscles in a relaxed state while keeping the tongue away from the bonding area. Designed with comfort and convenience in mind, these blocks are soft enough to trim with scissors if needed and feature openings to accommodate the Saliva Ejector (optional) and its companion Adapter Tube (optional). (20 per pack)</p> <ul style="list-style-type: none"> • Small • Medium | <p>300-032 300-033</p> |
|  | <p>Disposable Saliva Ejector and Adapter Tube Attachments for Disposable Tongue Block</p> <p>The Saliva Ejector is a perforated tube that's used to suction saliva from the patient's mouth. Constructed of medical grade vinyl tubing, it features a soft metal insert that enables the practitioner to shape it to best fit the individual patient. The companion Adapter Tube links the Saliva Ejector to the practitioner's vacuum system by a quick connect tip which speeds device changes. (3 per pack)</p> <ul style="list-style-type: none"> • Saliva Ejector • Adapter Tube | <p>300-034A 300-034B</p> |
|  | <p>Fuji ORTHO Gel Conditioner</p> <p>This gel conditioner is 20% polyacrylic acid, and is designed to mildly condition enamel. For easy application, the conditioner is available in syringes, with 25 disposable brush tips (3 syringes, 1.4 grams each, 25 tips).</p> | <p>480-116</p> |
|  | <p>Power Bond™ O.L.C. Light-Cured Bonding Conditioner</p> <p>Increases bond strength dramatically. Ideal for bonding to acrylic, composite or metal surfaces or any enamel surface where maximum adhesion is desired. Pre-mixed, one-step conditioner should be brushed on dry etched enamel prior to applying adhesive. This conditioner ensures a strong bond.</p> <p><i>Power Bond O.L.C. (7 grams)</i></p> | <p>471-062</p> |
|  | <p>Advantage Porcelain Bonding Conditioner</p> <p>Prepares dental porcelain for metal or ceramic bond brackets. Primer is a saline coupling agent which helps to chemically bond the adhesive to the porcelain.</p> <p>Porcelain Conditioner (16 cc)</p> | <p>480-013</p> |
|  | <p>Plastic Bracket Conditioner</p> <p>Prepares plastic brackets for direct bonding. Brush on base prior to applying adhesive. Provides a strong bond between the bracket base and adhesive. Must be refrigerated when not in use. (10 cc)</p> | <p>470-170</p> |