










# EXTRAORAL | J



ITEM	DESCRIPTION	CATALOG NO.																																												
	<p><b>Extraoral Laser Welded Facebow (Pro-Bow™)</b>                      This facebow is comfortable for the patient, plus it is easy to adjust and durable. The laser welded joint is much smoother and cleaner than a solder joint. The bow is made of solid wire with no hollow wire inserted. Hygienic, corrosion resistant, and more visually attractive. Permanently marked with a centerline. Wire diameter: .045; Angulation: 0°. Vertical loops on inner bow are adjustable.</p> <table border="1" data-bbox="427 646 1097 926"> <thead> <tr> <th>INNER ARCH</th> <th>SIZE NUMBER</th> <th>MM LENGTH</th> <th>CATALOG NO.</th> </tr> </thead> <tbody> <tr> <td>.045</td> <td>1</td> <td>83</td> <td>425-011</td> </tr> <tr> <td>.045</td> <td>2</td> <td>90</td> <td>425-012</td> </tr> <tr> <td>.045</td> <td>3</td> <td>97</td> <td>425-013</td> </tr> <tr> <td>.045</td> <td>4</td> <td>104</td> <td>425-014</td> </tr> <tr> <td>.045</td> <td>5</td> <td>111</td> <td>425-015</td> </tr> <tr> <td>.045</td> <td>6</td> <td>Straight</td> <td>425-016</td> </tr> </tbody> </table>	INNER ARCH	SIZE NUMBER	MM LENGTH	CATALOG NO.	.045	1	83	425-011	.045	2	90	425-012	.045	3	97	425-013	.045	4	104	425-014	.045	5	111	425-015	.045	6	Straight	425-016																	
INNER ARCH	SIZE NUMBER	MM LENGTH	CATALOG NO.																																											
.045	1	83	425-011																																											
.045	2	90	425-012																																											
.045	3	97	425-013																																											
.045	4	104	425-014																																											
.045	5	111	425-015																																											
.045	6	Straight	425-016																																											
	<p><b>Safety Neck Pads</b>                      These modules greatly reduce the risk of slingshot rebound from headgear. Reassembly after release is easy. The stainless steel coil spring allows the doctor to apply specific forces depending on which hole he inserts the facebow into.</p> <p>Neck Pad &amp; Safety Release with Stainless Steel Spring</p> <ul style="list-style-type: none"> <li>• SS Black Safety Neck Pad (Black Clip)</li> <li>• SS Black Safety Neck Pad (Glow-in-the-Dark Clip)</li> <li>• Black Pad Component (5 per pack)</li> </ul>	<p>426-002                      426-003                      426-000</p>																																												
	<p><b>Velpro™ Safety Neck Pad</b>                      Soft <i>Velpro</i> Material offers maximum patient comfort, while greatly reducing slingshot rebound from headgear. Reassembly after release is easy. The stainless steel coil spring allows the doctor to apply specific forces depending on which hole he or she inserts the facebow into. Comes as a complete assembly only with black clip.</p> <p>Black                      Purple                      Pink                      Blue                      Teal                      Black Clip Component (10 per pack)</p> <p><b>Force of Clip Component</b></p> <table border="1" data-bbox="427 1581 1097 1942"> <thead> <tr> <th>MARK NUMBER</th> <th>AVG. GRAMS</th> <th>AVG. OUNCES</th> <th>AVG. POUNDS</th> </tr> </thead> <tbody> <tr> <td>1</td> <td>291</td> <td>10</td> <td>0.65</td> </tr> <tr> <td>2</td> <td>414</td> <td>15</td> <td>0.98</td> </tr> <tr> <td>3</td> <td>524</td> <td>19</td> <td>1.17</td> </tr> <tr> <td>4</td> <td>613</td> <td>22</td> <td>1.37</td> </tr> <tr> <td>5</td> <td>730</td> <td>26</td> <td>1.63</td> </tr> <tr> <td>6</td> <td>875</td> <td>31</td> <td>1.95</td> </tr> <tr> <td>7</td> <td>966</td> <td>35</td> <td>2.16</td> </tr> <tr> <td>8</td> <td>1104</td> <td>39</td> <td>2.47</td> </tr> <tr> <td>9</td> <td>1196</td> <td>43</td> <td>2.67</td> </tr> <tr> <td>10</td> <td>1342</td> <td>49</td> <td>3.04</td> </tr> </tbody> </table>	MARK NUMBER	AVG. GRAMS	AVG. OUNCES	AVG. POUNDS	1	291	10	0.65	2	414	15	0.98	3	524	19	1.17	4	613	22	1.37	5	730	26	1.63	6	875	31	1.95	7	966	35	2.16	8	1104	39	2.47	9	1196	43	2.67	10	1342	49	3.04	<p>426-070                      426-071                      426-072                      426-073                      426-074                      426-020</p>
MARK NUMBER	AVG. GRAMS	AVG. OUNCES	AVG. POUNDS																																											
1	291	10	0.65																																											
2	414	15	0.98																																											
3	524	19	1.17																																											
4	613	22	1.37																																											
5	730	26	1.63																																											
6	875	31	1.95																																											
7	966	35	2.16																																											
8	1104	39	2.47																																											
9	1196	43	2.67																																											
10	1342	49	3.04																																											

ITEM	DESCRIPTION	CATALOG NO.
	<p><b>Face Mask (Reverse Pull)</b> For use in Class III cases, made of stainless steel and durable plastic. Features a lightweight, comfortable, design. Plus, comes in a convenient carrying case. Easily adjusts to fit all patients.</p>  <p>Blue Red Teal Pearl</p>	<p>426-055 426-056 426-057 426-058</p>
	<p><b>Extraoral Elastics</b> Available only in amber, ½" thickness (1000 per pack) 8 oz 16 oz</p>	<p>408-073 408-074</p>
	<p><b>Safety High-Pull Head Caps</b> Extra heavy duty construction for sturdy dependability. Easy to adjust.</p> <p><b>SS Black High Pull Cap</b> • Black Clip</p> <p><b>Safety High-Pull Head Cap Components</b> • Black Caps (5 per pack) • Black Clips (10 per pack)</p>	<p>426-012 426-010 426-020</p>
	<p><b>Class II Bow™</b> Designed for use with Class II elastics to distally drive or maintain anchorage of upper molars. Patient acceptance is much higher with <i>Class II Bow</i> over conventional headgear. Pre-soldered with pliable brass hooks for elastic use. Adjusting loops may be formed in wire, if desired. .045" wire diameter. (10 per pack)</p>	<p>100-090</p>
	<p><b>Chin Cup</b> Designed for use with either high-pull or vertical-pull headgear. Felt lined for comfort. Rigid plastic cup will not squeeze patient's chin or cheeks.</p>	<p>426-050</p>