

Bravo PrimoTM

Prescriptions and Ordering Information



MASEL[®]

Roth Prescription* (10 Per Pack)

BRACKET	DESCRIPTION	COLOR-CODE	TORQUE	ANG.	ROT.	M/D WIDTH	HGT.	O/G	SIDE	.018 ITEM NO.	.022 ITEM NO.	
MAXILLARY	Centrals	—	● Black	+12°	+5°	0°	0.130"	0.073"	0.122"	R	4499-315-10	4499-311-10
										L	4499-316-10	4499-312-10
	Laterals	—	● Pink	+8°	+9°	0°	0.109"	0.074"	0.122"	R	4499-325-10	4499-321-10
										L	4499-326-10	4499-322-10
	Cuspids	With Hook	● Green	-2°	+13°	4°M	0.117"	0.059"	0.122"	R	4499-335-10	4499-331-10
										L	4499-336-10	4499-332-10
		Without Hook		R	4499-358-10	4499-356-10						
				L	4499-359-10	4499-357-10						
	Bicuspids	With Hook	● Purple	-7°	0°	2°D	0.120"	0.065"	0.122"	R	4499-345-10	4499-341-10
										L	4499-344-10	4499-340-10
		Without Hook		R	4499-355-10	4499-351-10						
				L	4499-354-10	4499-350-10						
Bicuspids (Gingivly Offset Pads)	With Hook	● Purple	-7°	0°	2°D	0.120"	0.065"	0.122"	R	4499-745-10	4499-741-10	
									L	4499-744-10	4499-740-10	
	Without Hook		R	4499-755-10	4499-751-10							
			L	4499-754-10	4499-750-10							
MANDIBULAR	Anteriors	—	● Yellow	-1°	0°	0°	0.096"	0.077"	0.122"	R & L	4499-305-10	4499-301-10
	Cuspids	With Hook	● Blue	-11°	+7°	2°M	0.119"	0.068"	0.122"	R	4499-365-10	4499-361-10
										L	4499-366-10	4499-362-10
		Without Hook		R	4499-367-10	4499-363-10						
				L	4499-368-10	4499-364-10						
	Cuspids (OT)	With Hook	● Blue	0°	+3°	0°	0.119"	0.068"	0.122"	R	4499-561-10	
										L	4499-562-10	
		Without Hook		R	4499-563-10							
				L	4499-564-10							
	1st Bicuspids	With Hook	○ White	-17°	0°	4°D	0.119"	0.074"	0.122"	R	4499-377-10	4499-373-10
										L	4499-378-10	4499-374-10
		Without Hook		R	4499-375-10	4499-371-10						
				L	4499-376-10	4499-372-10						
	1st Bicuspids (Gingivly Offset Pads)	With Hook	○ White	-17°	0°	4°D	0.119"	0.074"	0.122"	R	4499-777-10	4499-773-10
										L	4499-778-10	4499-774-10
		Without Hook		R	4499-775-10	4499-771-10						
L				4499-776-10	4499-772-10							
2nd Bicuspids	With Hook	● Red	-22°	0°	4°D	0.119"	0.074"	0.122"	R	4499-387-10	4499-383-10	
									L	4499-388-10	4499-384-10	
	Without Hook		R	4499-385-10	4499-381-10							
			L	4499-386-10	4499-382-10							
2nd Bicuspids (Gingivly Offset Pads)	With Hook	● Red	-22°	0°	4°D	0.119"	0.074"	0.122"	R	4499-787-10	4499-783-10	
									L	4499-788-10	4499-784-10	
	Without Hook		R	4499-785-10	4499-781-10							
			L	4499-786-10	4499-782-10							

*MASEL'S® prescription does not imply endorsement by the doctor(s).

McLaughlin, Bennett, & Trevisi Prescription* (10 Per Pack)

BRACKET	DESCRIPTION	COLOR-CODE	TORQUE	ANG.	ROT.	M/D WIDTH	HGT.	O/G	SIDE	.022 ITEM NO.	
MAXILLARY	Centrals	—	 Black	+17°	+4°	0°	0.130"	0.073"	0.122"	R	4499-411-10
										L	4499-412-10
	Laterals	—	 Pink	+10°	+8°	0°	0.109"	0.074"	0.122"	R	4499-421-10
										L	4499-422-10
	Cuspids	With Hook	 Green	-7°	+8°	0°	0.117"	0.059"	0.122"	R	4499-431-10
										L	4499-432-10
		Without Hook		-7°	+8°	0°	0.117"	0.059"	0.122"	R	4499-456-10
										L	4499-457-10
	Cuspids (OT)	With Hook	 Green	0°	+8°	0°	0.117"	0.059"	0.122"	R	4499-531-10
										L	4499-532-10
		Without Hook		0°	+8°	0°	0.117"	0.059"	0.122"	R	4499-556-10
										L	4499-557-10
	Bicuspid	With Hook	 Purple	-7°	0°	0°	0.120"	0.065"	0.122"	R	4499-441-10
										L	4499-440-10
Without Hook		-7°		0°	0°	0.120"	0.065"	0.122"	R	4499-451-10	
									L	4499-450-10	
Bicuspid (Gingivly Offset Pads)	With Hook	 Purple	-7°	0°	0°	0.120"	0.065"	0.122"	R	4499-641-10	
									L	4499-640-10	
	Without Hook		-7°	0°	0°	0.120"	0.065"	0.122"	R	4499-651-10	
									L	4499-650-10	
MANDIBULAR	Anteriors	—	 Yellow	-6°	0°	0°	0.096"	0.077"	0.122"	R & L	4499-401-10
	Cuspids	With Hook	 Blue	-6°	+3°	0°	0.119"	0.068"	0.122"	R	4499-461-10
										L	4499-462-10
	Without Hook	-6°		+3°	0°	0.119"	0.068"	0.122"	R	4499-463-10	
									L	4499-464-10	
	Cuspids (OT)	With Hook	 Blue	0°	+3°	0°	0.119"	0.068"	0.122"	R	4499-561-10
										L	4499-562-10
		Without Hook		0°	+3°	0°	0.119"	0.068"	0.122"	R	4499-563-10
										L	4499-564-10
	1st Bicuspid	With Hook	 Gray	-12°	+2°	0°	0.119"	0.074"	0.122"	R	4499-473-10
										L	4499-474-10
		Without Hook		-12°	+2°	0°	0.119"	0.074"	0.122"	R	4499-471-10
										L	4499-472-10
	1st Bicuspid (Gingivly Offset Pads)	With Hook	 Gray	-12°	+2°	0°	0.119"	0.074"	0.122"	R	4499-673-10
L										4499-674-10	
Without Hook		-12°		+2°	0°	0.119"	0.074"	0.122"	R	4499-671-10	
									L	4499-672-10	
2nd Bicuspid	With Hook	 White	-17°	+2°	0°	0.119"	0.074"	0.122"	R	4499-483-10	
									L	4499-484-10	
	Without Hook		-17°	+2°	0°	0.119"	0.074"	0.122"	R	4499-481-10	
									L	4499-482-10	
2nd Bicuspid (Gingivly Offset Pads)	With Hook	 White	-17°	+2°	0°	0.119"	0.074"	0.122"	R	4499-683-10	
									L	4499-684-10	
	Without Hook		-17°	+2°	0°	0.119"	0.074"	0.122"	R	4499-681-10	
									L	4499-682-10	

*MASEL'S prescription does not imply endorsement by the doctor(s).

Roth Prescription* Single Cases $\frac{5-5}{5-5}$

BRACKETS PER CASE	DESCRIPTION	.018 ITEM NO.	.022 ITEM NO.
20	With Hooks on 3's, 4's & 5's	4499-018Q	4499-022Q
20	With Hooks on 3's, 4's & 5's and Gingivly Offset Pads on Lower 4's & 5's	4499-0180Q	4499-0220Q
20	With Hooks on 3's	4499-218Q	4499-222Q
20	With Hooks on 3's and Gingivly Offset Pads on Lower 4's & 5's	4499-2180Q	4499-2220Q
20	Without Hooks	4499-118Q	4499-122Q
20	Without Hooks and Gingivly Offset Pads on Lower 4's & 5's	4499-1180Q	4499-1220Q

McLaughlin, Bennett, & Trevisi Prescription* Single Cases $\frac{5-5}{5-5}$

BRACKETS PER CASE	DESCRIPTION	.022 ITEM NO.
20	With Hooks on 3's, 4's & 5's	4499-922Q
20	With Hooks on 3's, 4's & 5's and Gingivly Offset Pads on Lower 4's & 5's	4499-9220Q
20	With Hooks on 3's	4499-912Q
20	With Hooks on 3's and Gingivly Offset Pads on Lower 4's & 5's	4499-9120Q
20	Without Hooks	4499-902Q
20	Without Hooks and Gingivly Offset Pads on Lower 4's & 5's	4499-9020Q
20	With Hooks on 3's, 4's & 5's (0 Torque 3's)	4499-522Q
20	With Hooks on 3's, 4's & 5's (0 Torque 3's) and Gingivly Offset Pads on Lower 4's & 5's	4499-5220Q
20	With Hooks on 3's (0 Torque 3's)	4499-512Q
20	With Hooks on 3's (0 Torque 3's) and Gingivly Offset Pads on Lower 4's & 5's	4499-5120Q
20	Without Hooks (0 Torque 3's)	4499-502Q
20	Without Hooks (0 Torque 3's) and Gingivly Offset Pads on Lower 4's & 5's	4499-5020Q

*MASEL'S prescription does not imply endorsement by the doctor(s).



For more information, please contact us:

MaselOrtho.com

In the U.S.: 888.851.0533

Outside the U.S.: +(1) 760 448 8600

Canada: CERUM 800 661 9567

To email an order: Sales@Masel.net

Ortho Organizers, Inc.
1822 Aston Avenue
Carlsbad, CA 92008-7306
U.S.A.