

INTRODUCING  
**A REMARKABLE BREAKTHROUGH**  
IN CLASS III CORRECTION



**INTRODUCING THE ALL-NEW**  
**CARRIERE® MOTION™ CLASS III APPLIANCE**

The *Carriere Motion* Class III Appliance provides a new, remarkably easy-to-use and patient-friendly solution for Class III treatment. This discreet, comfortable appliance is direct bonded in just minutes, and is as easily tolerated as elastics alone! Imagine an appliance that gives you and your Class III patients an option without surgery or cumbersome, uncomfortable, and unsightly extra-oral devices. If you've ever struggled while tackling Class III cases, take a look at the all-new *Motion* Class III Appliance today!

Learn more about the *Motion* Class III Appliance at  
**888.851.0533** or **HenryScheinOrtho.com**.

Pre-Treatment



In Treatment



Post-Treatment



**HENRY SCHEIN®**  
ORTHODONTICS

The New Movement In Orthodontics™

# CARRIERE MOTION CLASS III APPLIANCE

The *Motion* Class III Appliance is a minimally invasive appliance intended to treat Class III malocclusions without extractions, painful orthognathic surgery, or traditionally uncomfortable and obtrusive face-masks.\* The *Motion* Class III Appliance is as easy to place as a bracket or buccal tube, and is designed to provide predictable results and create a harmonious, balanced profile, all while preserving the patient's natural features. The patient-friendly appliance is used at the beginning of treatment, prior to braces or aligners, and delivers rapid Anterior-Posterior correction to help reduce overall treatment times by several months. Intra-oral elastics connect the appliance with maxillary anchorage to activate the mandibular posterior segment and to move it bodily as a block for an ideal Class I malocclusion. Once the ideal platform is achieved, brackets or aligners are then placed to finish the case.

*Designed to provide predictable results and create a harmonious, balanced profile.*



KITS		ITEM NUMBER	
Motion Class III Standard Kit (1 set of each 23, 25, 27 mm)		424-407C	
DESCRIPTION/SIZE	COLOR-CODE	ITEM NUMBERS 1 pair	ITEM NUMBERS SINGLE (UNIVERSAL L/R)
Motion Class III 16 mm	White	424-416C	424-416
Motion Class III 18 mm	Blue	424-418C	424-418
Motion Class III 20 mm	Pink	424-420C	424-420
Motion Class III 23 mm	Yellow	424-423C	424-423
Motion Class III 25 mm	Red	424-425C	424-425
Motion Class III 27 mm	Green	424-427C	424-427
ACCESSORIES			
Typodont Motion Class III with Clear Aligner	-	-	631-033DNE
Typodont Motion Class III Maloccluded with Clear Aligner	-	-	631-034DNE
Motion Class III 5X Demonstration model	-	-	631-433X
Motion Elastics - Force 1, 1/4", 6 oz (50 packs of 100 elastics)	-	-	424-9F1
Motion Elastics - Force 2, 3/16", 8 oz (50 packs of 100 elastics)	-	-	424-9F2



The New Movement In Orthodontics™

For more information on our products and educational offerings, please contact us:  
**HenryScheinOrtho.com • Customer Service: 888.851.0533 • In the U.S.: 800.547.2000**  
**Outside the U.S.: +(1) 760 448 8600 • Canada: CERUM 800.661.9567**

Ortho Organizers, Inc. 1822 Aston Ave. Carlsbad, CA 92008-7306 USA

