


Maestro® Bracket Identification Chart



MAXILLARY

U1:	BLACK	DISTAL	GINGIVAL TIE WING
U2:	PINK	DISTAL	GINGIVAL TIE WING
U3:	GREEN	DISTAL	GINGIVAL TIE WING
U4 & 5:	PURPLE	DISTAL	GINGIVAL TIE WING
U6 & 7:	SEE ABOVE	TOP INDENT	

MANDIBULAR

L1 & 2:	YELLOW	GINGIVAL TIE WINGS = UNIVERSAL	
L3:	BLUE	GINGIVAL TIE WING	
L4 & 5:	SEE ABOVE		
L6 & 7:	SEE ABOVE	TOP INDENT	

ROTH

U6'S:	PURPLE
U7'S:	PINK
L4'S:	WHITE
L5'S:	RED
L6'S:	RED
L7'S:	BLUE

MCLAUGHLIN, BENNETT, AND TREVISI

U6'S:	PURPLE
U7'S:	GRAY
L4'S:	GRAY
L5'S:	WHITE
L6'S:	WHITE
L7'S:	YELLOW

ANDREWS

U6'S:	YELLOW
L6'S:	BLACK

LOW TORQUE

L6'S:	YELLOW
-------	--------