


Silhouette™ Poly-Sapphire™ Ceramic Brackets

MASEL®



Silhouette Poly-Sapphire Ceramic Brackets

The *Silhouette Poly-Sapphire* Ceramic Bracket gives you the clarity of a “sapphire” (monocrystalline) bracket, with less glare. The *Silhouette Poly-Sapphire* Ceramic Bracket provides the strength you require, while at the same time serving your patients’ needs by being resistant to stains and discoloration.


The *Silhouette Poly-Sapphire* Ceramic Bracket is the right choice for your practice.

- *Poly-Sapphire*, super-aesthetic appearance
 - An invisible look when you compare to sapphire brackets
 - The natural color of the tooth penetrates through the clear bracket
 - The glare commonly found with sapphire brackets is non-existent
- Hard surface
 - 35% less friction compared to sapphire brackets allows for more efficient tooth movement
 - Abrasion and high shock resistance
 - Resistant to stains and discoloration
- Low profile for patient comfort
- Mechanical bonding base for strong predictable bond strength and easy de-bonding
- Downdraft on tie-wings for easy, and safe ligation
- Single hook for a less bulky feel
- Available in both .018 and .022 slot sizes
- Available in both Roth and McLaughlin, Bennett, and Trevisi prescriptions
- Color-coded with an ID mark that can easily be removed



Also available for use as marketing tools to patients are the following typodonts and 15 times models:

Item Number	Description
2031-2886	Typodont - <i>Silhouette</i> (5 x 5)
2031-2886X	15X Demo Model - <i>Silhouette</i>



MAXILLARY

U1: BLACK - DISTAL-GINGIVAL TIE WING
 U2: PINK - DISTAL-GINGIVAL TIE WING
 U3: GREEN - DISTAL-GINGIVAL TIE WING
 U4 & 5: PURPLE - DISTAL-GINGIVAL TIE WING

MANDIBULAR

L1 & 2: YELLOW - GINGIVAL TIE WINGS = UNIVERSAL
 L3: BLUE - DISTAL-GINGIVAL TIE WING
 L4 & 5: SEE ABOVE

ROTH

LOWER 4'S: WHITE
 LOWER 5'S: RED

MCLAUGHLIN, BENNETT, AND TREVISI

LOWER 4'S: GRAY
 LOWER 5'S: WHITE

Standard Bracket Wheels Patient Packs


If you like the flexibility of loose brackets packaged in wheels, we can package them however you like. Each wheel can hold up to 10 cases, saving you space by holding larger quantities, or providing convenience for those who prefer to use one wheel per patient.

- Can hold 1 to 10 cases in each wheel



Less Friction Than Sapphire Brackets

Poly-Sapphire brackets experience 35% less friction than conventional sapphire brackets. Less friction optimizes the movement of teeth and shortens the overall patient treatment time as compared to regular sapphire brackets.



Why the Silhouette Bracket?



Mechanical bonding base
for strong predictable bond strength
and easy de-bonding



Low-Profile
for patient comfort



Single hook
for less bulky feel

Orthocosmetic Elastinol Preformed Archwires

This tooth-colored nickel titanium archwire has a special aesthetic coating that blends exceptionally well with ceramic or plastic brackets. **It won't stain or discolor, and resists cracking or chipping.** It facilitates fast tooth movement, is easy to engage even with the hardest-to-ligate bracket, and exerts a continuous, gentle force over a long period of time. This archwire is less brittle than other coated nickel titanium wires, enabling you to place a permanent stop in the posterior part of the wire.

Features and Benefits of Orthocosmetic Elastinol Preformed Archwires

- Longest lasting tooth-colored coating
- Super-elastic shape memory
- Recommended for initial treatment phases
- Round sizes only
- Ink center mark
- 5 per pack



Sizes*		Natural	
INCHES	MM	UPPER	LOWER
.014	.36	4833-114	4833-014
.016	.40	4833-116	4833-016
.018	.46	4833-118	4833-018

*Diameter includes coating

Orthocosmetic Preformed Ligature Ties

Perfect for Ceramic Brackets

A chemically inert, stain resistant tooth-colored coating gives these ties an "invisible" quality in the mouth. Your patients, especially those with ceramic or plastic brackets, will appreciate the improved aesthetics. Annealed dead-soft stainless steel gives the Orthocosmetic Preformed Ligature Ties just the right temper for fast, effortless ligating.

- 100 per pack

INCHES	MM	ITEM NO.
.012	.30	4910-122




Orthocosmetic Quik-Tie Ligatures

Aesthetic Coating Blends Naturally With Teeth—Perfect For Ceramic Brackets

MASEL's Orthocosmetic Quik-Tie™ Ligatures and Kobayashi Quik-Ties are available with a tooth-colored coating, and are easy to use. Just slip the wire over the bracket, twist to tighten, cut if needed, and tuck.

- White color
- Stain resistant
- Pre-cut, pre-tied to save you time
- 100 per pack



Orthocosmetic Quik-Tie

INCHES	MM	ITEM NO.
.010	.25	4911-210

Orthocosmetic Kobayashi Quik-Tie

INCHES	MM	ITEM NO.
.014	.36	4911-214

SILHOUETTE CERAMIC BRACKETS

Roth Prescription* Individual Brackets

BRACKET	DESCRIPTION	COLOR-CODE	TORQUE	ANG.	ROT.	M/D WIDTH	HGT.	O/G	SIDE	.018" ITEM NO.	.022" ITEM NO.	
MAXILLARY	Centrals	Black ●	+12°	+5°	0°	0.125"	0.075"	0.154"	R	2886-315	2886-311	
									L	2886-316	2886-312	
	Laterals	Pink ●	+8°	+9°	0°	0.125"	0.089"	0.152"	R	2886-325	2886-321	
									L	2886-326	2886-322	
	Cuspids	with Hook	Green ●	0°	+11°	+4°M	0.125"	0.075"	0.202"	R	2886-335	2886-331
										L	2886-336	2886-332
		without Hook								R	2886-358	2886-356
										L	2886-359	2886-357
	Bicuspid	with Hook	Purple ●	-7°	0°	0°	0.125"	0.075"	0.202"	R	2886-345	2886-341
										L	2886-344	2886-340
without Hook		R								2886-355	2886-351	
		L								2886-354	2886-350	
MANDIBULAR	Anteriors	Yellow ●	0°	0°	0°	0.12"	0.091"	0.151"	R & L	2886-305	2886-301	
	Cuspids	with Hook	Blue ●	-11°	+5°	+2°M	0.125"	0.075"	0.202"	R	2886-365	2886-361
										L	2886-366	2886-362
		without Hook								R	2886-367	2886-363
										L	2886-368	2886-364
	1st Bicuspid	with Hook	White ○	-17°	0°	0°	0.125"	0.076"	0.202"	R	2886-377	2886-373
										L	2886-378	2886-374
		without Hook								R	2886-375	2886-371
										L	2886-376	2886-372
	2nd Bicuspid	with Hook	Red ●	-22°	0°	0°	0.125"	0.078"	0.202"	R	2886-387	2886-383
L										2886-388	2886-384	
without Hook		R								2886-385	2886-381	
		L								2886-386	2886-382	

Roth Prescription* Case Set-ups in Wheels

SET-UP	BRACKETS PER CASE	DESCRIPTION	.018" ITEM NO.	.022" ITEM NO.
<u>5-5</u> 5-5	20	with Hooks on 3's, 4's, & 5's	2886-0461	2886-0471
<u>5-5</u> 5-5	20	with Hooks on 3's	2886-0481	2886-0491
<u>5-5</u> 5-5	20	without Hooks	2886-0441	2886-0451
<u>5-5</u> 5-5	20 (10 Cases)	with Hooks on 3's, 4's, & 5's	2886-0461-10	2886-0471-10
<u>5-5</u> 5-5	20 (10 Cases)	with Hooks on 3's	2886-0481-10	2886-0491-10
<u>5-5</u> 5-5	20 (10 Cases)	without Hooks	2886-0441-10	2886-0451-10
<u>5-5</u>	10	with Hooks on 3's, 4's, & 5's	2886-0462	2886-0472
<u>5-5</u>	10	with Hooks on 3's	2886-0482	2886-0492
<u>5-5</u>	10	without Hooks	2886-0442	2886-0452
<u>3-3</u> <u>3-3</u>	12	with Hooks on 3's	2886-036	2886-037
<u>3-3</u> <u>3-3</u>	12	without Hooks	2886-034	2886-035
<u>3-3</u>	6	with Hooks on 3's	2886-032	2886-033
<u>3-3</u>	6	without Hooks	2886-030	2886-031

SILHOUETTE CERAMIC BRACKETS

McLaughlin, Bennett, and Trevisi Prescription* Individual Brackets

BRACKET	DESCRIPTION	COLOR-CODE	TORQUE	ANG.	ROT.	M/D WIDTH	HGT.	O/G	SIDE	.018" ITEM NO.	.022" ITEM NO.	
MAXILLARY	Centrals	Black ●	+17°	+4°	0°	0.125"	0.075"	0.154"	R	2886-415	2886-411	
									L	2886-416	2886-412	
	Laterals	Pink ●	+10°	+8°	0°	0.125"	0.089"	0.152"	R	2886-425	2886-421	
									L	2886-426	2886-422	
	Cuspids	with Hook	Green ●	-7°	+8°	0°	0.125"	0.075"	0.202"	R	2886-435	2886-431
										L	2886-436	2886-432
		without Hook		R	2886-458	2886-456						
				L	2886-459	2886-457						
	Bicuspid	with Hook	Purple ●	-7°	0°	0°	0.125"	0.075"	0.202"	R	2886-445	2886-441
										L	2886-444	2886-440
without Hook		R		2886-455	2886-451							
		L		2886-454	2886-450							
MANDIBULAR	Anteriors	Yellow ●	-6°	0°	0°	0.12"	0.091"	0.151"	R & L	2886-405	2886-401	
	Cuspids	with Hook	Blue ●	-6°	+3°	0°	0.125"	0.075"	0.202"	R	2886-465	2886-461
										L	2886-466	2886-462
		without Hook		R	2886-467	2886-463						
				L	2886-468	2886-464						
	1st Bicuspid	with Hook	Gray ●	-12°	0°	0°	0.125"	0.076"	0.202"	R	2886-477	2886-473
										L	2886-478	2886-474
		without Hook		R	2886-475	2886-471						
				L	2886-476	2886-472						
	2nd Bicuspid	with Hook	White ●	-17°	0°	0°	0.125"	0.078"	0.202"	R	2886-487	2886-483
L										2886-488	2886-484	
without Hook		R		2886-485	2886-481							
		L		2886-486	2886-482							

McLaughlin, Bennett, and Trevisi Prescription* Case Set-ups in Wheels

SET-UP	BRACKETS PER CASE	DESCRIPTION	.018" ITEM NO.	.022" ITEM NO.
$\frac{5-5}{5-5}$	20	with Hooks on 3's, 4's, & 5's	2886-9461	2886-9471
$\frac{5-5}{5-5}$	20	with Hooks on 3's	2886-9481	2886-9491
$\frac{5-5}{5-5}$	20	without Hooks	2886-9441	2886-9451
$\frac{5-5}{5-5}$	20 (10 Cases)	with Hooks on 3's, 4's, & 5's	2886-9461-10	2886-9471-10
$\frac{5-5}{5-5}$	20 (10 Cases)	with Hooks on 3's	2886-9481-10	2886-9491-10
$\frac{5-5}{5-5}$	20 (10 Cases)	without Hooks	2886-9441-10	2886-9451-10
$\frac{5-5}{5-5}$	10	with Hooks on 3's, 4's, & 5's	2886-9462	2886-9472
$\frac{5-5}{5-5}$	10	with Hooks on 3's	2886-9482	2886-9492
$\frac{5-5}{5-5}$	10	without Hooks	2886-9442	2886-9452
$\frac{3-3}{3-3}$	12	with Hooks on 3's	2886-936	2886-937
$\frac{3-3}{3-3}$	12	without Hooks	2886-934	2886-935
$\frac{3-3}{3-3}$	6	with Hooks on 3's	2886-932	2886-933
$\frac{3-3}{3-3}$	6	without Hooks	2886-930	2886-931



MASEL[®]

Ortho Organizers, Inc.
1822 Aston Avenue
Carlsbad, CA 92008
U.S.A.
MaselOrtho.com

To Order by Phone

Domestic

Ph 888 851 0533 760 448 8600
Fax 800 888 7244 760 448 8602
Email Sales@Masel.net

International

Ph 760 448 8600
Fax 760 448 8607
Email Sales@Masel.net

Canada

Ph +(1) 800 661 9567
Fax +(1) 800 361 5088
Email cerumsup@telusplanet.net

Australia

Ph +1300 ORDERZ or +612 9540 6111
Fax (02) 9524 7623 or +612 9524 7623
Email sales@naol.com.au