

Masque Ceramic Brackets

Looking for an aesthetic bracket that resists breakage, and is easy to debond? *MASEL*® is pleased to offer you the *Masque* Ceramic Bracket, a translucent bracket made from polycrystalline alumina material. In addition, it blends perfectly with the natural tooth tone to give your patients the aesthetic appearance they desire.

Choose *Masque* Ceramic Brackets and you'll enjoy these features and benefits:

- High purity, polycrystalline alumina material for maximum strength and high flexibility to prevent breakage
- Translucent material blends perfectly with natural tooth tone for an aesthetically pleasing look
- Mechanical bonding base, strong predictable bond strength, and easy de-bonding
- Generous tie-wing area for safe and easy ligation
- Universal hooks on 4's and 5's reduce inventory requirements
- Available in both Roth and McLaughlin, Bennett, and Trevisi prescriptions
- Color-coded with an ID mark that can easily be removed



MASQUE CERAMIC BRACKETS

Roth Prescription* Individual Brackets

BRACKET	DESCRIPTION	COLOR-CODE	TORQUE	ANG.	ROT.	M/D WIDTH	HGT.	O/G	SIDE	.018" ITEM NO.	.022" ITEM NO.	
MAXILLARY	Centrals	Black ●	+12°	+5°	0°	0.125"	0.075"	0.154"	R	2887-315	2887-311	
									L	2887-316	2887-312	
	Laterals	Pink ●	+8°	+9°	0°	0.125"	0.089"	0.152"	R	2887-325	2887-321	
									L	2887-326	2887-322	
	Cuspids	with Hook	Green ●	0°	+11°	+4°M	0.125"	0.075"	0.202"	R	2887-335	2887-331
										L	2887-336	2887-332
		without Hook								R	2887-358	2887-356
										L	2887-359	2887-357
	Bicuspids	with Hook	Purple ●●	-7°	0°	0°	0.125"	0.075"	0.202"	R	2887-345	2887-341
										L	2887-344	2887-340
without Hook		R								2887-355	2887-351	
		L								2887-354	2887-350	
MANDIBULAR	Anteriors	Yellow ●	0°	0°	0°	0.12"	0.091"	0.151"	R & L	2887-305	2887-301	
	Cuspids	with Hook	Blue ●	-11°	+5°	+2°M	0.125"	0.075"	0.202"	R	2887-365	2887-361
										L	2887-366	2887-362
		without Hook								R	2887-367	2887-363
										L	2887-368	2887-364
	1st Bicuspids	with Hook	White ○○	-17°	0°	0°	0.125"	0.076"	0.202"	R	2887-377	2887-373
										L	2887-378	2887-374
		without Hook								R	2887-375	2887-371
										L	2887-376	2887-372
	2nd Bicuspids	with Hook	Red ●●	-22°	0°	0°	0.125"	0.078"	0.202"	R	2887-387	2887-383
L										2887-388	2887-384	
without Hook		R								2887-385	2887-381	
		L								2887-386	2887-382	

Roth Prescription* Case Set-ups in Wheels

SET-UP	BRACKETS PER CASE	DESCRIPTION	.018" ITEM NO.	.022" ITEM NO.
<u>5-5</u> 5-5	20	with Hooks on 3's, 4's, & 5's	2887-0461	2887-0471
<u>5-5</u> 5-5	20	with Hooks on 3's	2887-0481	2887-0491
<u>5-5</u> 5-5	20	without Hooks	2887-0441	2887-0451
<u>5-5</u> 5-5	20 (10 Cases)	with Hooks on 3's, 4's, & 5's	2887-0461-10	2887-0471-10
<u>5-5</u> 5-5	20 (10 Cases)	with Hooks on 3's	2887-0481-10	2887-0491-10
<u>5-5</u> 5-5	20 (10 Cases)	without Hooks	2887-0441-10	2887-0451-10
<u>5-5</u>	10	with Hooks on 3's, 4's, & 5's	2887-0462	2887-0472
<u>5-5</u>	10	with Hooks on 3's	2887-0482	2887-0492
<u>5-5</u>	10	without Hooks	2887-0442	2887-0452
<u>3-3</u> <u>3-3</u>	12	with Hooks on 3's	2887-036	2887-037
<u>3-3</u> <u>3-3</u>	12	without Hooks	2887-034	2887-035
<u>3-3</u>	6	with Hooks on 3's	2887-032	2887-033
<u>3-3</u>	6	without Hooks	2887-030	2887-031

* MASEL's prescription does not imply endorsement by the doctors.

MASQUE CERAMIC BRACKETS

McLaughlin, Bennett, and Trevisi Prescription* Individual Brackets

BRACKET	DESCRIPTION	COLOR-CODE	TORQUE	ANG.	ROT.	M/D WIDTH	HGT.	O/G	SIDE	.018" ITEM NO.	.022" ITEM NO.	
MAXILLARY	Centrals	Black ●	+17°	+4°	0°	0.125"	0.075"	0.154"	R	2887-415	2887-411	
									L	2887-416	2887-412	
	Laterals	Pink ●	+10°	+8°	0°	0.125"	0.089"	0.152"	R	2887-425	2887-421	
									L	2887-426	2887-422	
	Cuspids	with Hook	Green ●	-7°	+8°	0°	0.125"	0.075"	0.202"	R	2887-435	2887-431
										L	2887-436	2887-432
		without Hook		R	2887-458	2887-456						
				L	2887-459	2887-457						
	Bicuspids	with Hook	Purple ●●	-7°	0°	0°	0.125"	0.075"	0.202"	R	2887-445	2887-441
										L	2887-444	2887-440
without Hook		R		2887-455	2887-451							
		L		2887-454	2887-450							
MANDIBULAR	Anteriors	Yellow ●	-6°	0°	0°	0.12"	0.091"	0.151"	R & L	2887-405	2887-401	
	Cuspids	with Hook	Blue ●	-6°	+3°	0°	0.125"	0.075"	0.202"	R	2887-465	2887-461
										L	2887-466	2887-462
		without Hook		R	2887-467	2887-463						
				L	2887-468	2887-464						
	1st Bicuspids	with Hook	Gray ●●	-12°	0°	0°	0.125"	0.076"	0.202"	R	2887-477	2887-473
										L	2887-478	2887-474
		without Hook		R	2887-475	2887-471						
				L	2887-476	2887-472						
	2nd Bicuspids	with Hook	White ●●●	-17°	0°	0°	0.125"	0.078"	0.202"	R	2887-487	2887-483
L										2887-488	2887-484	
without Hook		R		2887-485	2887-481							
		L		2887-486	2887-482							

McLaughlin, Bennett, and Trevisi Prescription* Case Set-ups in Wheels

SET-UP	BRACKETS PER CASE	DESCRIPTION	.018" ITEM NO.	.022" ITEM NO.
$\frac{5-5}{5-5}$	20	with Hooks on 3's, 4's, & 5's	2887-9461	2887-9471
$\frac{5-5}{5-5}$	20	with Hooks on 3's	2887-9481	2887-9491
$\frac{5-5}{5-5}$	20	without Hooks	2887-9441	2887-9451
$\frac{5-5}{5-5}$	20 (10 Cases)	with Hooks on 3's, 4's, & 5's	2887-9461-10	2887-9471-10
$\frac{5-5}{5-5}$	20 (10 Cases)	with Hooks on 3's	2887-9481-10	2887-9491-10
$\frac{5-5}{5-5}$	20 (10 Cases)	without Hooks	2887-9441-10	2887-9451-10
$\frac{5-5}{5-5}$	10	with Hooks on 3's, 4's, & 5's	2887-9462	2887-9472
$\frac{5-5}{5-5}$	10	with Hooks on 3's	2887-9482	2887-9492
$\frac{5-5}{5-5}$	10	without Hooks	2887-9442	2887-9452
$\frac{3-3}{3-3}$	12	with Hooks on 3's	2887-936	2887-937
$\frac{3-3}{3-3}$	12	without Hooks	2887-934	2887-935
$\frac{3-3}{3-3}$	6	with Hooks on 3's	2887-932	2887-933
$\frac{3-3}{3-3}$	6	without Hooks	2887-930	2887-931

* MASEL's prescription does not imply endorsement by the doctors. © 2012 Ortho Organizers, Inc. All rights reserved.