















Glossary Symbols

Symbols	Title & Designation No. of Standard	Title of Symbols	Reference No	Explanatory Text
	ISO 15223-1	Do not re-use	5.4.2	Indicates a medical device that is intended for one use, or for use on a single patient during a single procedure
	ISO 15223-1	Caution	5.4.4	Indicates the need for the user to consult the instruction for use for important cautionary information such as warnings and precautions that cannot, for a variety of reasons, be presented on the medical device itself.
	ISO 15223-1	Consult instructions for use	5.4.3	Indicates the need for the user to consult the instructions for use.
	ISO 15223-1	Authorized representative in the European Community	5.1.2	Indicates the Authorized Representative in the European Community.
	ISO 15223-1	Manufacturer	5.1.1	Indicates the medical device manufacturer, as defined in EU Directives 90/385/EEC, 93/42/EEC and 98/79/EEC.
	ISO 15223-1	Batch code	5.1.5	Indicates the manufacturer's batch code so that the batch or lot can be identified.
	ISO 15223-1	Use-by date	5.1.4	Indicates the date after which the medical device is not to be used.
	ISO 15223-1	Temperature limit	5.3.7	Indicates the temperature limits to which the medical device can be safely exposed.
	ISO 15223-1	Keep away from sunlight	5.3.2	Indicates a medical device that needs protection from light sources.
	ISO 15223-1	Catalog number	5.1.6	Indicates the manufacturer's catalog number so that the medical device can be identified.
	ISO 15223-1	Contains or presence of natural rubber latex	5.4.5	Indicates the presence of natural rubber or dry natural rubber latex as a material of construction within the medical device or the packaging of a medical device
Rx Only	21 CFR 801.109(b)(1)	Prescription device labeling statement	N/A	Caution: Federal law restricts this device to sale by or on the order of a (licensed healthcare practitioner)
	GHS Hazard Pictogram	Corrosive	N/A	Indicates a corrosive substance that will destroy and damage other substances with which it comes into contact.

	GHS Hazard Pictogram	Flammable	N/A	Indicates a flammable substance that can be ignited very easily and will burn very vigorously when ignited.
	GHS Hazard Pictogram	Irritant	N/A	Indicates a material that can cause immediate health effects such as skin rashes or respiratory irritation.